


SELECTION GUIDE FOR FUME HOODS AND VENTILATION SYSTEMS FOR SCIENCE EDUCATION



Protecting your
laboratory environment

LABCONCO

THE PURPOSE OF A FUME HOOD:

To capture, contain and remove fumes and contaminants generated by scientific experiments.

VENTILATION REFERENCES AND STANDARDS

Hood Location:

According to the current ANSI Z9.5 recommendations, fume hoods shall not be placed in a location that will adversely affect their performance. Placement with respect to open windows, doorways, personal traffic flow areas and supply air diffusers can influence containment. A good rule of thumb is at least 10 feet from a door or window and at least 5 feet from a supply air diffuser.

ASHRAE 110-1995

American Society of Heating, Refrigerating and Air-Conditioning Engineers

- Three part test: Smoke generation, face velocity profile, tracer gas release @ 4 liters/minute.
- Rated As Manufactured (AM), As Installed (AI) and As Used (AU).

Industrial Ventilation-ACGIH

American Conference of Governmental Industrial Hygienists

- Stack discharge height 1.3-2.0 times the building height

NFPA 45

National Fire Protection Association

- Materials of construction should have flame spread of 25 or less
- Fume hood exhaust air should not be recirculated

ASTM E84-01

American Society for Testing and Materials

- Test method to determine flame spread
 - Mineral fiber cement board flame spread 0
 - Red oak flooring flame spread 100

ANSI Z9.5-2003

- All hoods shall be equipped with a flow indicator, flow alarm or face velocity monitor
- Ductless hoods should only be used with non-hazardous material

OSHA 29 CFR Part 1910

Occupational Health and Safety

- 2.5 linear feet of hood space per person

DEMONSTRATION HOODS

Advanced

Protector® ClassMate® Hood



Requires Ductwork
and Blower

*Protector ClassMate
Laboratory Hood
6970401 shown with
SpillStopper Work Surface
6986401 and Protector
Acid Storage Cabinet
9901000*

- For general chemistry applications in an instructional setting
- Clear tempered safety glass back and sides and taller front viewing window provide enhanced visibility
- Ducts to the outside

All Protector ClassMate models listed will be equipped with 4 service fixtures (gas, air, vacuum and water) and 2 electrical duplex receptacles.

Size (w x d x h)

- ___ 48" x 32.25" x 59" Labconco #6970401
- ___ 60" x 32.25" x 59" Labconco #6970501
- ___ 72" x 32.25" x 59" Labconco #6970601

ACCESSORIES:

SpillStopper™ Work Surface (6" x 3" cupsink cutout left rear)

- ___ 4' Labconco #6986401
- ___ 5' Labconco #6986501
- ___ 6' Labconco #6986601

6" x 3" Cup Sink

- ___ Labconco #4005200

Airflow Monitor, 115VAC

- ___ Guardian Jr. (Analog), Labconco #9743204
- ___ Guardian 335 Digital, Labconco #9743203

Blower

- ___ Fiberglass
- ___ Coated Steel

See Blower Selection Guide on page 8. Consult Labconco for help with blower sizing. Ductwork length and configuration are needed to size the blower correctly.

For base cabinet or stand selection, see page 6.

Additional models are available. Contact Labconco for details and ordering information.

DEMONSTRATION HOODS

Intermediate

Protector® XVS Ventilation Station



Requires Ductwork and Blower or FilterMate Portable Exhauster

Protector XVS Ventilation Station 4863010 shown with Work Surface 3908401 and Base Stand 3746711

- For light general chemistry in educational applications
- Clear tempered safety glass sides and top, steel back
- Ducts to the outside. Can recirculate into the room with FilterMate™ Portable Exhauster

Size (w x d x h)

- ___ 24" x 22.81" x 22.75" Labconco #4862010
- ___ 36" x 22.81" x 22.75" Labconco #4863010
- ___ 48" x 22.81" x 22.75" Labconco #4864010

ACCESSORIES:

Solid Epoxy Dished Work Surface

- ___ 2' Labconco #3908400
- ___ 3' Labconco #3908401
- ___ 4' Labconco #3908402

Airflow Monitor, 115VAC

- ___ Guardian 500 (Analog), Labconco #3944700
- ___ Guardian 1000 (Digital), Labconco #3908800

Fluorescent Light Kit 115VAC

- ___ Labconco #3909200

Blower or Filtered Exhauster

- ___ Remote Coated Steel Blower, Labconco #4863500
- ___ FilterMate Portable Exhauster, Labconco #3970001

Organic Filter for FilterMate Portable Exhauster

- ___ Labconco #3923400

See www.labconco.com/services/request-chemical-assessment to complete a Chemical Assessment form for your application. Other filter types are also available.

Duct Transition Adapters

- ___ Upper to 6" diameter duct, Labconco #3912401
- ___ Lower to 6" diameter duct, Labconco #3912403

For stand or base cabinet selection, see page 6.

Additional models are available. Contact Labconco for details and ordering information.

DEMONSTRATION HOODS

Simple

Protector® Demonstration Hood



Requires Ductwork and Blower or FilterMate Portable Exhauster

Protector Demonstration Hood System 3945000

- For very light chemistry, chemical reactions, classroom observation
- Tempered safety glass sides, top and back
- Ducts to the outside. Can recirculate into room with FilterMate™ Portable Exhauster

Size (w x h x d)

- ___ 36" x 22.75" x 29.1" Labconco #3944300
 - ___ 36" x 32" x 29.1" Labconco #3944301
 - ___ Protector Demonstration Hood System (22.75"h) Labconco #3945000
 - ___ Protector Demonstration Hood System (32"h) Labconco #3945001
- Protector Demonstration Hood System includes Enclosure, Epoxy Work Surface, Cart and FilterMate Portable Exhauster.

ACCESSORIES:

Solid Epoxy Dished Work Surface

- ___ 3' Labconco #3908401

Airflow Monitor, 115VAC

- ___ Guardian 500 (Analog), Labconco #3944700
- ___ Guardian 1000 (Digital), Labconco #3908800

Fluorescent Light Kit 115VAC

- ___ Labconco #3909200

Blower or Filtered Exhauster

- ___ Remote Coated Steel Blower, Labconco #4863500
- ___ FilterMate Portable Exhauster, Labconco #3970001

Organic Filter for FilterMate Portable Exhauster

- ___ Labconco #3923400

See www.labconco.com/services/request-chemical-assessment to complete a Chemical Assessment form for your application. Other filter types are also available.

Duct Transition Adapters

- ___ Upper to 6" diameter duct, Labconco #3912401
- ___ Lower to 6" diameter duct, Labconco #3912403

For stand or base cabinet selection, see page 6.

EXPERIMENTAL HOODS

Advanced

Protector® Premier® Hood



Requires Ductwork
and Blower

4' Protector Premier Laboratory Hood 100400002 is shown with SpillStopper Work Surface 9500410 and Protector Standard Storage Cabinet 9900000.

- For general chemistry applications
- One-piece molded fiberglass liner is corrosion resistant, durable, fire resistant and light reflective

All Protector Premier models listed will be equipped with 2 service fixtures and 1 electrical duplex receptacle

Size (w x d x h)

- ___ 48.0" x 31.7" x 59.0" Labconco #100400002
- ___ 60.0" x 31.7" x 59.0" Labconco #100500002
- ___ 72.0" x 31.7" x 59.0" Labconco #100600002
- ___ 96.0" x 31.7" x 59.0" Labconco #100800002

ACCESSORIES:

SpillStopper™ Work Surface (6" x 3" sink cut out left rear)

- ___ 4' Labconco #9500410
- ___ 5' Labconco #9500510
- ___ 6' Labconco #9500610
- ___ 8' Labconco #9500810

6" x 3" Cup Sink

- ___ Labconco #4005200

Cold Water Gooseneck

- ___ Labconco #9827900

Airflow Monitor, 115VAC

- ___ Guardian, Labconco #9413300
- ___ Guardian Digital, Labconco #9413400

Blower

- ___ Fiberglass
- ___ Coated Steel

See Blower Selection Guide on page 8. Consult Labconco for help with blower sizing. Ductwork length and configuration are needed to size the blower correctly.

For base cabinet or stand selection, see page 6.

Additional models are available. Contact Labconco for details and ordering information.

EXPERIMENTAL HOODS

Intermediate

Fiberglass 30 Hood



Requires Ductwork

Fiberglass 30 Laboratory Hood 3030000 shown with Work Surface 4882806 and Protector Standard Storage Cabinet 9900200

- For light to moderate chemistry applications
- One-piece molded fiberglass liner is corrosion resistant, durable, fire resistant and light reflective

The Fiberglass 30 Hood model listed is equipped with a built-in exhaust blower. Electrical receptacles and service fixtures are sold as accessories.

Size (w x d x h)

- ___ 30" x 30.2" x 48" Labconco #3030000

ACCESSORIES:

Solid Epoxy Flat Work Surface (no sink cut out)

- ___ 30" Labconco #4882806

Electrical Duplex Kit, 115VAC, 20amp

- ___ Labconco #2834800

Service Fixture Kit, Gas

- ___ Labconco #2832605

Airflow Monitor, 115VAC

- ___ Guardian 500 (Analog), Labconco #9743213
- ___ Guardian (Digital), Labconco #9743201

See Blower Selection Guide on page 8. Consult Labconco for help with blower sizing. Ductwork length and configuration are needed to size the blower correctly.

For base cabinet or stand selection, see page 6.

Additional models are available. Contact Labconco for details and ordering information.

EXPERIMENTAL HOODS

Simple

Basic 47 and 70 Hoods



Requires Ductwork

Basic 47 Laboratory Hood 2247300 shown with Work Surface 4882803 and Protector Standard Storage Cabinet 9900000

- For non-corrosive chemistry applications
- Epoxy-coated steel interior and exterior
- Basic Hoods are less deep and high than Premier and ClassMate models. Designed to fit existing counter tops and cramped areas.

All Basic Hood models listed are equipped with a built-in exhaust blower. Blower can be relocated to the roof. Electrical receptacles and service fixtures are sold as accessories.

Size (w x d x h)

- ___ 47" x 25" x 53" Labconco #2247300
- ___ 70" x 25" x 53" Labconco #2246300

ACCESSORIES:

Solid Epoxy Flat Work Surface (no sink cut out)

- ___ 4' Labconco #4882803
- ___ 6' Labconco #4882804

6" x 3" Cup Sink & Turret Kit (mounts on right side wall of hood)

- ___ Labconco #4852500

Electrical Duplex Kit, 115VAC, 20amp

- ___ Labconco #4844400

Service Fixture Kit, Gas

- ___ Labconco #2832605

Airflow Monitor, 115VAC

- ___ Guardian Jr. (Analog), Labconco #9743208
- ___ Guardian (Digital), Labconco #9743207

For base cabinet or stand selection, see page 6.

Additional models are available. Contact Labconco for details and ordering information.

DUAL ACCESS HOODS

Protector® Pass-Through Hood



Requires Ductwork and Blower

Protector Pass-Through Laboratory Hood 9942502 shown with Work Surface 9972500 and Protector Standard Storage Cabinets 9900200

- Vertical rising sashes front and back permit dual access
- Audible alarm prevents simultaneous opening of both sashes
- Chemical resistant composite panel liner
- Available with right side viewing window for classrooms or without viewing window for through-wall installations

All models include one electrical duplex receptacle.

Size (w x d x h)

With right side viewing window:

- ___ 48" x 33.2" x 59" Labconco #9942402
- ___ 60" x 33.2" x 59" Labconco #9942502
- ___ 72" x 33.2" x 59" Labconco #9942602

Without viewing window:

- ___ 48" x 33.2" x 59" Labconco #9942403
- ___ 60" x 33.2" x 59" Labconco #9942503
- ___ 72" x 33.2" x 59" Labconco #9942603

ACCESSORIES:

SpillStopper™ Work Surface (no sink cut out)

- ___ 4' Labconco #9972400
- ___ 5' Labconco #9972500
- ___ 6' Labconco #9972600

Airflow Monitor, 115VAC

- ___ Guardian 500 (Analog), Labconco #9743213
- ___ Guardian 1000 (Digital), Labconco #9743211

Blower

- ___ Fiberglass
- ___ Coated Steel

See Blower Selection Guide on page 8. Consult Labconco for help with blower sizing. Ductwork length and configuration are needed to size the blower correctly.

For base cabinet or stand selection, see page 6.

DUCTLESS HOODS

Paramount® Ductless Enclosure



Paramount Ductless Enclosure 6963302 shown with Work Surface 3909901 and Telescoping Base Stand with casters 3746711

- Requires no ducting and can be moved from lab to lab
- For light to moderate chemistry applications
- Clear tempered safety glass sides, back and baffle
- Uses carbon filters to trap or treat chemical fumes
- Uses HEPA filters to trap particulates
- Equipped with a built-in blower

Size (w x d x h)

- ___ 36" x 29" x 53" Labconco #6963302
- ___ 48" x 29" x 53" Labconco #6963402

ACCESSORIES:

Paramount Filters

- ___ Ammonia-Amine Filter Labconco #6938202
- ___ Acid-Sulfur Filter Labconco #6938201
- ___ Organic Vapor Filter Labconco #6938200
- ___ Formaldehyde Filter Labconco #6938203
- ___ Mixed Bed Filter (four chemicals listed above) Labconco#6938204
- ___ HEPA Filter Labconco #6938100
- See www.labconco.com/services/request-chemical-assessment to complete a Chemical Assessment form for your application.

Solid Epoxy Dished Work Surface

- ___ 3' Labconco #3909901
- ___ 4' Labconco #3909902

Hydraulic Lift Base Stands

Adjust from 25.5" to 33.5" for easy access to filters and mobility through 80" doorways. Stands include leveling feet which can be replaced with 5" diameter toe-locking casters.

- ___ 3' Electric Labconco #3780300
- ___ 3' Manual Labconco #3780400
- ___ 4' Electric Labconco #3780301
- ___ 4' Manual Labconco #3780401
- ___ Casters for field replacement of leveling feet Labconco #3784000

For additional stand or base cabinet selection, see column to right. Additional models are available. Contact Labconco for details and ordering information.

BASE CABINETS



Standard Storage Cabinet 9900000

Standard Storage Cabinet

- ___ 24" single right hinge Labconco #9900300
- ___ 24" single left hinge Labconco #9900600
- ___ 30" dual doors Labconco #9900200
- ___ 36" dual doors Labconco #9900100
- ___ 48" dual doors Labconco #9900000

Acid Storage Cabinet

- ___ 24" single right hinge Labconco #9901300
- ___ 24" single left hinge Labconco #9901500
- ___ 30" dual doors Labconco #9901200
- ___ 36" dual doors Labconco #9901100
- ___ 48" dual doors Labconco #9901000
- ___ Vent Kit for Acid Storage Cabinet Labconco #3591100

Solvent Storage Cabinet

- ___ 24" single right hinge Labconco #9902300
- ___ 24" single left hinge Labconco #9902400
- ___ 30" dual doors Labconco #9902200
- ___ 36" dual doors Labconco #9902100
- ___ 48" dual doors Labconco #9902000
- ___ Vent Kit for Solvent Storage Cabinet Labconco #9910000

BASE STANDS



Telescoping Base Stand 3746702


Telescoping Base Stands

Adjust to seven height positions from 27.5" to 33.5" in 1" increments

- ___ 2' Labconco #3767000 (Use with Protector XVS)
- ___ 2' with casters Labconco #3767010 (Use with Protector XVS)
- ___ 3' Labconco #3746701 (Use with Paramount Ductless Enclosure, Protector XVS, Demonstration Hood and Fiberglass 30)
- ___ 3' with casters Labconco #3746711 (Use with Paramount Ductless Enclosure, Protector XVS and Demonstration Hood)
- ___ 4' Labconco #3746702 (Use with Paramount Ductless Enclosure, Protector XVS, Protector Hoods and Basic 47)
- ___ 4' with casters Labconco #3746712 (Use with Paramount Ductless Enclosure and Protector XVS)
- ___ 5' Labconco #3746703 (Use with Protector Hoods)
- ___ 6' Labconco #3746704 (Use with Protector Hoods and Basic 70)

TYPICAL VENTILATION SYSTEM

This diagram details the many components that are needed to complete a typical fume hood ventilation system. All of these components, as well as others, are available from Labconco.



DUCTWORK AND ACCESSORIES

Thermoplastic Duct, 10 feet lengths

- ___ 6" Labconco #4708600
- ___ 8" Labconco #4718900
- ___ 10" Labconco #7027200
- ___ 12" Labconco #5602000



Duct Couplings, Female

- ___ 6" Labconco #4708900
- ___ 8" Labconco #4719200
- ___ 10" Labconco #7027500
- ___ 12" Labconco #5602300



Duct Couplings, Male

- ___ 6" Labconco #2144700
- ___ 8" Labconco #4719900
- ___ 10" Labconco #7027800
- ___ 12" Labconco #7067300



Elbows, 90 degree

- ___ 6" Labconco #4708700
- ___ 8" Labconco #4719000
- ___ 10" Labconco #7027300
- ___ 12" Labconco #5602100



Elbows, 45 degree

- ___ 6" Labconco #4708800
- ___ 8" Labconco #4719100
- ___ 10" Labconco #7027400
- ___ 12" Labconco #5602200



Duct Reducers

- ___ 6" x 8" Labconco #5605900
- ___ 8" x 10" Labconco #5606000
- ___ 10" x 12" Labconco #5606100



Zero Pressure Weather Caps

- ___ 6" Labconco #4722200
- ___ 8" Labconco #4722300
- ___ 10" Labconco #7095100
- ___ 12" Labconco #5622100



Blower Transitions, for Blowers 7068000-7068700

- ___ 8" Labconco #4722400
- ___ 10" Labconco #4722401



Flexible Duct Connections

- ___ for 6" Duct, Labconco #4861800
- ___ for 8" Duct, Labconco #4726500
- ___ for 10" Duct, Labconco #7034200
- ___ for 12" Duct, Labconco #5621400



Manual Duct Dampers

- ___ 6" Labconco #4724200
- ___ 8" Labconco #4741300
- ___ 10" Labconco #5983400
- ___ 12" Labconco #5981200



BLOWERS

Labconco Blowers have been specifically designed to meet the demanding day-to-day ventilation requirements for a wide range of hood choices. Coated Steel Blowers are recommended for low to moderately corrosive applications. They are offered in a broad spectrum of capacities from 300 to 2150 CFM at static pressures to 1.50" H₂O.


Labconco Fiberglass Blowers are ideal for fume hood exhaust systems in moderate to highly corrosive conditions. They are offered in a broad range of capacities from 250 to 3500 CFM at static pressures to 4.0".

Coated Steel Blowers and Fiberglass Blowers feature forward curved impellers to ensure quiet operation and optimum air delivery. They are belt-driven with adjustable sheave and self-adjusting gravity belt tightener. The variable exhaust orientation adjusts to eight outlet positions for installation flexibility (vertical discharge is recommended).

How To Calculate Equivalent Length of Hood Duct

The equivalent length duct method is used to properly size blowers. To determine the correct blower, first determine the total equivalent length of duct in the installation.

| Elbow Equivalent Resistance Chart (In Feet Of Straight Duct) | | | | |
|---|---------|---------|----------|----------|
| | 6" Duct | 8" Duct | 10" Duct | 12" Duct |
| 90° Elbow | 12' | 15' | 20' | 25' |
| 45° Elbow | 6' | 7.5' | 10' | 12.5' |



Calculation:
 Straight Length,
 20' + 5' + 3'
 = 28'
 Weathercap
 = 5'
 90° Elbow
 = 25'
 Total Equivalent
 Feet
 = 58'
 (in straight duct)



Protecting your
laboratory environment

LABCONCO®

Labconco Corporation
 8811 Prospect Avenue
 Kansas City, MO 64132-2696
 816-333-8811 • 800-821-5525
 FAX: 816-363-0130
 Home Page: www.labconco.com
 E-Mail: labconco@labconco.com



Blower Type/Equivalent Duct Run

| Hood | Sash Open | Face Velocity (fpm) | Airflow (CFM) | Duct Run Dia.(in) | Steel 25' | Fiberglass 25' | Steel 50' | Fiberglass 50' | Steel 75' | Fiberglass 75' | Steel 100' | Fiberglass 100' |
|--|-----------|---------------------|---------------|-------------------|-----------|----------------|-----------|----------------|-----------|----------------|------------|-----------------|
| 4 ¹ Protector Premier Hood | 100% | 100 | 725 | 12 | 7068800 | 7181000 | 7068800 | 7181000 | 7068800 | 7181000 | 7068800 | 7181000 |
| 4 ¹ Protector Premier Hood | 60% | 100 | 435 | 12 | 7068000 | 7180000 | 7068000 | 7180000 | 7068000 | 7180000 | 7068000 | 7180000 |
| 5 ¹ Protector Premier Hood | 100% | 100 | 955 | 12 | 7069000 | 7181200 | 7069000 | 7181200 | 7069000 | 7181200 | 7069000 | 7181200 |
| 5 ¹ Protector Premier Hood | 60% | 100 | 575 | 12 | 7068800 | 7180800 | 7068800 | 7180800 | 7068800 | 7180800 | 7068800 | 7180800 |
| 6 ¹ Protector Premier Hood | 100% | 100 | 1180 | 12 | 7069000 | 7181400 | 7069000 | 7181400 | 7069200 | 7181600 | 7069200 | 7181600 |
| 6 ¹ Protector Premier Hood | 60% | 100 | 710 | 12 | 7068800 | 7181000 | 7068800 | 7181000 | 7068800 | 7181000 | 7068800 | 7181000 |
| 4 ¹ Protector ClassMate Hood | 60% | 100 | 570 | 8 | 7068200 | 7180200 | 7068200 | 7180200 | 7068200 | 7180400 | 7068400 | 7180400 |
| 5 ¹ Protector ClassMate Hood | 60% | 100 | 730 | 8 | 7068400 | 7180400 | 7068400 | 7180600 | 7068600 | 7180600 | 7068600 | 7180600 |
| 6 ¹ Protector ClassMate Hood | 60% | 100 | 890 | 10 | 7068600 | 7180600 | 7068600 | 7180600 | 7068600 | 7180600 | 7068600 | 7180600 |
| 4 ¹ Protector Pass-Through Hood | 100% | 100 | 830 | 10 | 7068400 | 7180400 | 7068400 | 7180400 | 7068600 | 7180600 | 7068600 | 7180600 |
| 4 ¹ Protector Pass-Through Hood | 60% | 100 | 530 | 8 | 7068000 | 7180200 | 7068200 | 7180200 | 7068200 | 7180200 | 7068200 | 7180200 |
| 5 ¹ Protector Pass-Through Hood | 100% | 100 | 1100 | 12 | 7069000 | 7181400 | 7069200 | 7181400 | 7069200 | 7181400 | 7069200 | 7181400 |
| 5 ¹ Protector Pass-Through Hood | 60% | 100 | 700 | 8 | 7068200 | 7180400 | 7068400 | 7180400 | 7068400 | 7180600 | 7068600 | 7180600 |
| 6 ¹ Protector Pass-Through Hood | 100% | 100 | 1350 | 12 | 7069200 | 7181400 | 7069400 | 7181400 | 7069400 | 7181600 | 7069400 | 7181600 |
| 6 ¹ Protector Pass-Through Hood | 60% | 100 | 865 | 10 | 7069000 | 7180600 | 7068600 | 7180600 | 7068600 | 7180600 | 7068600 | 7180600 |
| 2 ¹ Protector XVS Ventilation Station | 100% | 80 | 95 | 6 | 4863500 | — | 4863500 | — | 4863500 | — | 4863500 | — |
| 3 ¹ Protector XVS Ventilation Station | 100% | 80 | 145 | 6 | 4863500 | — | 4863500 | — | 4863500 | — | 4863500 | — |
| 4 ¹ Protector XVS Ventilation Station | 100% | 80 | 200 | 6 | 4863500 | — | 4863500 | — | 4863500 | — | 4863500 | — |
| Protector Demonstration Hood, 22.75" tall | 100% | 80 | 145 | 6 | 4863500 | — | 4863500 | — | 4863500 | — | 4863500 | — |
| Protector Demonstration Hood, 32" tall | 100% | 80 | 185 | 6 | 4863500 | — | 4863500 | — | 4863500 | — | 4863500 | — |