

pH, Indicators & Dyes



pH, Indicators & Dyes Kit

TEACHING STRATEGY

Kit activities are designed for introductory (grades 6-8) and intermediate-level / high school (grades 8-9) students. There are a total of 6 discrete Model investigations, 3 Inquiry Path investigations, and 2 Going Further investigations.

ACTIVITY 1 can be done independently with introductory-level students. ACTIVITIES 2, 3 and 4 are linked - each precedes the other and are recommended for advanced middle school or high school-level students.



ACTIVITY 1 pH, Indicators & Dyes


This kit supports 40 students working in groups of 4

KIT CONTENTS

- 40 medicine cups
- 10 cotton fabric (4" squares)
- 10 plastic jars , with lids (2oz)
- 100 commercial filter paper strips (1 x 4")
- 1 box, (72) microscope slides
- 1 pH paper strips, pkg./100 (0 to 14)
- 10 pipets
- 10 clothespins
- 1 sheet, black construction paper
- ball, twine (50')

CD-ROM

- ▶ TEACHER GUIDE set PDF
- ▶ STUDENT GUIDE set PDF
- ACTIVITY 1 *Learning About pH, Indicators, and Dyes* PDF
- ACTIVITY 2 *Learning About the Chemistry of Plant Pigments* PDF
- ACTIVITY 3 *Learning About Making pH Paper Test Strips* PDF
- ACTIVITY 4 *Learning About Testing and Accuracy* PDF
- ▶ BACKGROUND INFORMATION set PDF
 - pH, INDICATORS, & DYES (INTRODUCTORY) PDF
 - PIGMENT EXTRACTION (INTRODUCTORY) PDF
 - PLANT PIGMENT BIOMOLECULES & THEIR STRUCTURE (INTERMEDIATE) PDF
 - pH & POH (INTERMEDIATE) PDF
 - PLANT PIGMENTS AS pH INDICATORS (INTRODUCTORY) PDF
 - ALTERING PLANT PIGMENT BIO MOLECULES (INTERMEDIATE) PDF
 - LEARNING HOW LITMUS PAPER IS MADE (INTERMEDIATE) PDF
- ▶ *Learning About pH* PPT and .MOV (with *Study Notes* (PDF))



NEEDED BUT NOT SUPPLIED

From the grocery store

- Bag, Blueberries (bag, fresh or frozen) ACTIVITIES: 1,2
- Box, Baking soda (sodium bicarbonate) ACTIVITY 1
- Bottle, White vinegar (acetic acid) ACTIVITIES: 1, 2, 3
- 10+ Bottles water (8 oz; 240mL), or distilled water ACTIVITIES: 1,2,3
- Container, lemon / lime juice concentrate ACTIVITIES: 1,2, 3
- Container, carbonated water, apple juice ACTIVITIES: 1, 2, 3
- Container, table salt or rock salt ACTIVITIES: 1, 2, 3
- Container, liquid household detergent, hand soap ACTIVITIES: 1, 2, 3
- Container, liquid household ammonia / bleach ACTIVITIES: 1, 2, 3
- Package, Ammonium Aluminum Sulfate¹ - optional ACTIVITY: 1
- Package, Drano® ACTIVITY: 2

From a local source / school / home

- 10 Pairs of scissors ACTIVITIES: 1, 3
- 10 Plastic spoons ACTIVITIES: 1
- 20 Sheets, white copy paper ACTIVITIES: 1
- 10 Coffee filters ACTIVITIES: 1
- 10 Sheets, acid-free paper (fine stationary paper) ACTIVITIES: 1
- 10 Sets, colored pencils / colored markers ACTIVITIES: 1
- 10 Sheets, watercolor paper (optional) ACTIVITIES: 1
- 10 Microscope slides ACTIVITY 2b
- 10 Marking pens ACTIVITIES: 1, 2

ACCESSING DIGITAL CONTENT

- FILE TYPES:
- .pdf - portable digital file (Mac / PC)
 - notes.pdf - powerpoint pdf with lecture notes
 - .pps - powerpoint (PC)
 - .mov - Quick Time
 - iPod .mov Quick Time (iPod / iPad)

NOTE: Animations play within PDF files

See: *Content Delivery Information* to learn more!

INFORMATION TABLE 1		
ACTIVITY / Level	Model Investigations	Inquiry Path Investigations
1 INTRODUCTORY	1a <i>Blueberries as Biological Indicators</i> (40 minutes)	<i>Dye Color Design</i> (40 minutes)
	1b <i>Blueberries as a Fabric Colorant</i> (40 minutes)	<i>The Blueberry Cyanidin Molecule as a pH Indicator</i> (40 minutes)
		GOING FURTHER: <i>Exploring Natural Dyes in Fabrics</i> (40 minutes)
2 INTERMEDIATE	2a <i>The Blueberry Cyanidin Molecule as a pH Indicator</i>	<i>Determine pH range for a biopigment</i> (40 minutes)
	2b <i>Converting a Chlorophyll Molecule to a Phaeophytin Molecule</i> (30 minutes)	
3 INTERMEDIATE	<i>Making pH Paper Test Strips</i> (40 minutes)	GOING FURTHER: <i>Exploring Colored Paper as pH Indicators</i> (40 minutes)
4 INTERMEDIATE	<i>Evaluating pH Paper Test Strips</i> (40 minutes)	

PREP & EXPERIMENT TIMES

INFORMATION TABLE 2 Expected prep and Experimental Run Times		
ACTIVITY	PREP	LAB
1a <i>Blueberries as Biological Indicators</i>	10 minutes	40 minutes
1b <i>Blueberries as a Fabric Colorant</i>	20 minutes	40 minutes
INQUIRY PATH: <i>Dye Color Design</i>	20 minutes	40 minutes
INQUIRY PATH: <i>The Blueberry Cyanidin Molecule as a pH Indicator</i>	20 minutes	40 minutes
GOING FURTHER: <i>Exploring Natural Dyes in Fabrics</i>	20 minutes	40 minutes
2a <i>The Blueberry Cyanidin Molecule as a pH Indicator</i>	30 minutes	30 minutes
2b <i>Converting a Chlorophyll Molecule to a Phaeophytin Molecule</i>	20 minutes	30 minutes
INQUIRY PATH: <i>Determine pH range for a Bio Pigment</i>	20 minutes	40 minutes
3 <i>Making pH Paper Test Strips</i>	20 minutes	40 minutes
GOING FURTHER: <i>Exploring Colored Paper as pH Indicators</i>	20 minutes	40 minutes
<i>Evaluating pH Paper Test Strips</i>	10 minutes	40 minutes

ACTIVITY 1 INTRODUCTORY INVESTIGATION

Learning About pH, Indicators, and Dyes

OVERVIEW

In this structured investigation, students learn about biological pigments; how they can be used as pH indicators and be manipulated as fabric colorants.

OBJECTIVES

- ▶ *To understand* the difference between a biological pigment, a dye, and other kinds of pigments as chemical compounds.
- ▶ *To understand* what solvents are, and how biological pigments can be extracted.
- ▶ *To understand* the concept of pH and the role of biological indicators.
- ▶ *To design extraction methods for harvesting biological pigments in water-based solutions.*
- ▶ *To use extracted biological pigments as indicators.*
- ▶ *To dye fabric swatches.*

MATERIALS (Per Class)

In the Kit:

- 40 medicine cups
- 10 cotton fabric squares
- 10 plastic jars, with lids (2oz)
- 10 clothespins
- ball, twine (50')
- CD-ROM

Needed but not supplied:

From the grocery store

- 1 bag, Blueberries (bag, fresh or frozen)
- 1 box, Baking soda (sodium bicarbonate)
- 1 bottle, White vinegar (glacial acetic acid)
- 10+ bottles water (8 oz; 240mL) - or distilled water

Needed but not supplied (cont.)

- 1 container, lemon / lime juice concentrate (STRONG ACID examples)
- 1 container, carbonated water, apple juice (WEAK ACID examples)
- 1 container, table salt or rock salt (VERY WEAK ACID examples)
- 1 container, liquid household detergent, hand soap (WEAK BASE examples)
- 1 container, liquid household ammonia / bleach (STRONG BASE examples)
- 1 Bottle, Ammonium Aluminum Sulfate² - optional

From a local source / school / home:

- 10 pairs of scissors
- 10 plastic spoons
- 20 sheets, white copy paper
- 10 coffee filters
- 10 sets, colored pencils / colored markers

BEGINNINGS ...

▶ Have your students read the checked information topics:

Background Information Topics

- ✓ pH, INDICATORS, & DYES
- ✓ PIGMENT EXTRACTION?
- PLANT PIGMENT BIOMOLECULES & THEIR STRUCTURE
- pH & POH
- PLANT PIGMENTS AS pH INDICATORS
- ALTERING PLANT PIGMENT BIO MOLECULES
- LEARNING HOW LITMUS PAPER IS MADE

▶ Have students research the answers to these questions ...
(Students can obtain answers using *Background Information* and other web resources, as you direct.)

TEACHERS NOTE: These topics are available as individual PDF files on the CD-ROM. They can be emailed or uploaded to the school server. *Background* material is available to students as individual topic PDFs.

- ✓ Define and give an example of a universal solvent.
Water. Most chemical compounds can be dissolved in it.
- ✓ Is a pigment a dye? organize your answer in an information table.
In some ways, and not in others:

Information Table 1		
Characteristic	Dye	Pigment
Chemical Compound (consisting of 2 or more atoms)	✓	✓
Easily Dissolves (is soluble in ...)	✓	
Present in living things		✓ (biological pigment)
Appear colored due to light reflection	✓	✓
Man-made colored substance	✓	✓
Natural colored substance	✓	✓

- ✓ What is the bio pigment in blueberries?
The red/blue flavonoid bio pigment cyanidin.
- ✓ What is the usual purpose of biological pigment extraction?
To obtain a colored, transparent solution.
- ✓ What is an indicator solution?
A solution that indicates, through a color change, the presence of an acid or alkali (acid or base).
- ✓ What is pH; what are acids and alkalis?
pH is a numerical measure of how acidic or basic a solution and/or substance is. pH is measured on a pH scale: (0) STRONG ACID to (14) STRONG BASE.

An acid is a chemical compound that has a pH less than 7.

A base is a chemical compound having a pH greater than 7.

- ✓ Can extracted biological pigments be both dyes and indicators?
Yes.
- ✓ What is vinegar and baking soda? What are their chemical formulae?
vinegar is a 5-7% solution of acetic acid CH_3COOH .
baking soda is sodium bicarbonate NaHCO_3
- ✓ What is the pH of vinegar and baking soda?
vinegar is pH 3
baking soda is pH 9
- ✓ If an acid is reacted with a base, what happens?
A chemical compound called a salt is produced, along with water.

GETTING STARTED ...

► **MODEL INVESTIGATION 1a** - *Blueberries as biological indicators* introduces students to:

- *Extraction* of a biological pigment (cyanidin, an anthocyanin in the skin of blueberries - See: *Information Table 2*) to create a *colored, transparent, water solution*
- *Visual recognition* that this extracted pigment is an indicator - changes color with the change in pH
- Students will use this model investigation as a guide in designing an extraction method to produce a biological dye to color a fabric a color (at a specific pH).

► **MODEL INVESTIGATION 1b** - *Blueberries as a Fabric Colorant* introduces students to:

- *Dyeing a fabric swatch using an extracted biological pigment.*
- Students will use this model investigation as a guide in designing a method to dye fabric a specific color.

PREP NOTES

► Provide the following materials and work areas:

- ✓ bottle of white vinegar (weak acid)
- ✓ box of baking soda (weak base)
- ✓ set up a 'drying area' by installing a string 'clothesline' in the classroom for student use

► Provide the following materials for each student group:

- ✓ 4 medicine cups
- ✓ 3 blueberries (fresh or frozen)
- ✓ 1 plastic jar (2oz), filled with bottled or distilled water[▲]
- ✓ 1 bottle water (8oz; 240mL)
- ✓ 1 plastic spoon
- ✓ 1 sheet, copy paper
- ✓ 1 fabric swatch
- ✓ 1 scissors

[▲] Do not use tap water which has chlorine, it can interfere with the extraction

► **(INVESTIGATION 1a)** Following 15 minutes of blueberry anthocyanin extraction, visit each group and add:

- Cup #1 a pinch of baking soda (enough to change the pH to dark blue)
- Cup #2 nothing
- Cup #3 1 capful (3mL) vinegar (enough to change the pH to a pink color)

TEACHERS NOTE

One key action is that students observe the leaching of the bio pigment, cyanidin, directly from the skin of the blueberry - without any application of heat, or physical maceration. They should conclude that water is the solvent - forming a colored solution. These same blueberries are then set aside - ready for use in the INQUIRY PATH.