

Synergy™ H1 Hybrid Multi-Mode Microplate Reader

Synergy H1 is a flexible monochromator-based multi-mode microplate reader that can be turned into a high-performance hybrid system with the addition of a filter-based optical module. The monochromator optics uses a third generation quadruple grating design that allows working at any excitation or emission wavelength with a 1 nm step. This system supports top and bottom fluorescence intensity, UV-visible absorbance and high performance luminescence detection. It is the ideal system for all the standard microplate applications found in life science research laboratories. The filter module is a completely independent add-on that includes its own light source, and a high performance dichroic-based wavelength selection system. With its very high optical efficiency, this module supports advanced detection modes such as Fluorescence Polarization, Time-Resolved Fluorescence and filtered luminescence (e.g. BRET). A dual reagent injection system is available to automate inject/read assays such as ion channels assays or flash luminescence assays (e.g. luciferase or ATP assays).



Compatible with Take3™ Multi-Volume Plate with sixteen 2 μ L microspots for low-volume assays.



Features:

- Patent-Pending Hybrid Technology combines flexible monochromator detection with high performance dichroic-based filter detection
- Compatible with Take3™ Multi-Volume Plate: Samples down to 2 μ L volume can be measured. Especially useful when working with precious samples, for fast and accurate DNA/RNA quantification at 260 nm
- Quadruple grating monochromator for maximum flexibility and ease of use
- Dichroic-based filter optics, for best performance and advanced detection technologies such as fluorescence polarization and time-resolved fluorescence
- Comes with Gen5™ Data Analysis Software: reader control, advanced data analysis and flexible Excel export in one software package



Quadruple grating monochromator system: Ease of use and flexibility.



Easy-to-use filter system with magnetic filter cubes that can be swapped in a matter of seconds.

Models:

- Synergy H1m: Monochromator-based
- Synergy H1f: Filter-based
- Synergy H1mf: Hybrid

Dual reagent dispenser option available with all models.

Accessories:

- Take3™ Multi-Volume Plate
- BioStack™: 30 or 50 plate stacker
- Gen5™ Secure for 21CFR part 11 compliance
- Product Qualification Package
- Fluorescence, Luminescence and Absorbance Test Plates

Specifications:

Fluorescence Intensity

| | |
|-----------------------|--|
| Light Source: | Xenon Flash Lamp |
| Wavelength Range: | Monochromators: 250-700 nm (900 nm option) Filters: 200-700 nm (900 nm option) |
| Wavelength Selection: | Double grating monochromators (Top/Bottom) and deep blocking bandpass filters / dichroic mirrors (Top) |
| Sensitivity Top: | Monochromators: fluorescein 2.5 pM typical (0.25 fmol/well 384-well plate) Filters/mirrors: fluorescein 1 pM typical (0.1 fmol/well 384-well plate) |

Luminescence

| | |
|--------------------|---|
| Wavelength Range: | 300-700 nm |
| Dynamic Range: | > 6 Decades |
| Sensitivity (ATP): | Monochromator system: 20 amol ATP typical (flash) Filter system: 10 amol ATP typical (flash) |

Absorbance

| | |
|-----------------------|---------------------------------|
| Light Source: | Xenon Flash Lamp |
| Wavelength Range: | 230-999 nm, 1 nm increment |
| Wavelength Selection: | Monochromator |
| Bandwidth: | 4 nm (230-285 m), 8 nm (>285 m) |
| Measurement Range: | 0 to 4.0 OD |
| OD Accuracy: | < 1 % at 2 OD typical |
| OD Precision: | < 0.5 % at 2 OD typical |

Fluorescence Polarization

| | |
|-----------------------|---|
| Wavelength Range: | 300-700 nm (900 nm option) |
| Wavelength Selection: | Deep blocking bandpass filters/dichroic mirrors (top) |
| Sensitivity Top: | 5 mP at 1 nM fluorescein typical |

Time Resolved Fluorescence

| | |
|-----------------------|--|
| Light Source: | Xenon Flash Lamp |
| Wavelength Range: | 200-700 nm (900 nm option) |
| Wavelength Selection: | Deep blocking bandpass filters/dichroic mirrors |
| Sensitivity Top: | Europium 100 fM typical (10 amol/well in 384-well plate) |

Dispensers

| | |
|----------------------|-------------------------------|
| Number of Injectors: | 2 syringe pumps |
| Dispense Volume: | 5 - 1000 µL in 1 µL increment |
| Dead Volume: | 1 mL, 100 µL with backflush |

Speed (Minimum Kinetic Interval)

| | |
|-----------|------------|
| 96-well: | 11 seconds |
| 384-well: | 22 seconds |

General

| | |
|---|---|
| Microplate Types: | 1- to 384-well plates Compatible with Take3™ Multi-Volume Plate with 2 µL microspots |
| Shaking: | Yes |
| Dimensions: | 15.4"W x 18.6"D x 12.9"H (39.1 x 47.2 x 32.8 cm) |
| Weight: | 50 lbs (22.5 kg) |
| Temperature Control: | RT to 45°C |
| Top Optics Adjustment: | Automated |
| BioStack Compatible (Automation-Ready): | Yes |
| Software: | Gen5™ Data Analysis Software |

*Specifications subject to change



For customer service, call 1-800-766-7000.
To fax an order, use 1-800-926-1166.
To order online: www.fishersci.com
©2010 Thermo Fisher Scientific Inc. All rights reserved.