



Biological Buffers and Ultra Pure Reagents

One Call. One Source.
A World of Ultra Pure Biochemicals.



Are MP Buffers
in your corner?


Find it at fishersci.com

 **fisher scientific**
part of Thermo Fisher Scientific

Theoretical Considerations

Since buffers are essential for controlling the pH in many biological and biochemical reactions, it is important to have a basic understanding of how buffers control the hydrogen ion concentration. Although a lengthy, detailed discussion is impractical, some explanation of the buffering phenomena is important.

Let us begin with a discussion of the equilibrium constant (K) for weak acids and bases. Acids and bases which do not completely dissociate in solution, but instead exist as an equilibrium mixture of undissociated and dissociated species, are termed weak acids and bases. The most common example of a weak acid is acetic acid. In solution, acetic acid exists as an equilibrium mixture of acetate ions, hydrogen ions, and undissociated acetic acid. The equilibrium between these species may be expressed as follows:



where k_1 is the dissociation rate constant of acetic acid to acetate and hydrogen ions and k_2 is the association rate constant of the ion species to form acetic acid. The rate of dissociation of acetic acid, $-d(\text{HAc})/dt$, may be expressed by the following equation:

$$- \frac{d(\text{HAc})}{dt} = k_1 (\text{HAc})$$

which shows the rate of dissociation to be dependent upon the rate constant of dissociation (k_1) and the concentration of acetic acid (HAc).

Similarly, the rate association, $d(\text{HAc})/dt$, which is dependent upon the rate constant of association (k_2) and the concentration of acetate and hydrogen ions, may be shown as:

$$\frac{d(\text{HAc})}{dt} = k_2 (\text{H}^+) (\text{Ac}^-)$$

Since, under equilibrium conditions, the rates of dissociation and association must be equal, they may be expressed as:

$$k_1 (\text{HAc}) = k_2 (\text{H}^+) (\text{Ac}^-)$$

Or

$$\frac{k_1}{k_2} = \frac{(\text{H}^+) (\text{Ac}^-)}{(\text{HAc})}$$

If we now let $k_1/k_2 = K_a$, the equilibrium constant, the equilibrium expression becomes:

$$K_a = \frac{(\text{H}^+) (\text{Ac}^-)}{(\text{HAc})}$$

which may be rearranged to express the hydrogen ion concentration in terms of the equilibrium constant and the concentrations of undissociated acetic acid and acetate ions as follows:

$$(\text{H}^+) = K_a \frac{(\text{HAc})}{(\text{Ac}^-)}$$

Since pH is defined as $-\log(\text{H}^+)$, if the equilibrium expression is converted to $-\log$:

$$-\log(\text{H}^+) = -\log K_a - \log \frac{(\text{HAc})}{(\text{Ac}^-)}$$

And by substituting pH and pK_a :

$$\text{pH} = pK_a - \log \frac{(\text{HAc})}{(\text{Ac}^-)}$$

Or

$$\text{pH} = pK_a + \log \frac{(\text{Ac}^-)}{(\text{HAc})}$$

When the concentration of acetate ions equals the concentration of acetic acid, $\log (Ac^-)/(HAc)$ becomes zero, and the pH equals pK_a . As a result, the pK_a of a weak acid or base generally indicates the pH of the center of the buffering region.

pK_a values are generally determined by titration. The free acid of the material to be measured is carefully titrated with a suitable base, and using a calibrated automatic recording titrator, the titration curve is recorded. A general titration curve for a typical monobasic weak acid is shown in Figure 1. The point of inflection indicates the pK_a value.




Figure 1.
Typical Titration Curve of a weak acid

Using acetic acid as an example, it has now been demonstrated that $pH = pK_a$ when the concentrations of acetic acid and acetate ions are equal. This buffering action helps explain how the hydrogen ion concentration (H^+) remains relatively unaffected by external influences. Let's look at a hypothetical buffer system, HA ($pK_a = 7.000$) and (A^-). If we consider a non-buffered system to which a strong acid is added, we can observe a significant change in pH. For example, if 100 mL of 1.000×10^{-2} M HCl are added to 1.000 liter of 1.000 M NaCl at pH 7.000, the hydrogen ion concentration (H^+)_f of the final 1.100 liters of solution may be calculated by:

$$(H^+)_f \times Vol_f = (H^+)_i \times Vol_i$$

$$(H^+)_f \times 1.100 = 1.000 \times 10^{-7} \times 0.100$$

$$(H^+)_f = 9.09 \times 10^{-9}$$

$$-\log (H^+)_f = -\log (9.09 \times 10^{-9})$$

$$pH = 3.04$$

Thus, it can be observed that the addition of 1.0×10^{-3} moles of hydrogen ion to the unbuffered system resulted in a change in pH from 7.000 to 3.04.

Now, using the hypothetical buffer system, a 1.000 M solution of HA at pH 7.000 can be shown initially as:

$$(HA) = (A^-) = 0.500 \text{ M}$$

$$pH = pK + \log \frac{(A^-)}{(HA)}$$

$$pH = 7.000 + \log \frac{0.500}{0.500}$$

$$pH = 7.000$$

If we add to this system 100 mL of 1.000×10^{-2} M HCl, 1.000×10^{-3} moles of A must be converted to 1.000×10^{-3} moles of HA.

The resulting equation thus becomes:

$$pH = 7.000 + \log \frac{0.499/1.100}{0.501/1.100}$$

$$pH = 7.000 - 0.002$$

$$pH = 6.998$$

So it can be seen that in the buffered system the pH has changed by only 0.002 pH units, compared to a change of almost 4 pH units in the unbuffered system.

In summary, the principles involved in hydrogen ion buffer systems have been very basically illustrated. Beginning with an understanding of equilibrium, pH and pK_a , we have attempted to demonstrate how buffering capacity is determined and how a buffered system may effectively resist changes in pH.

Practical Considerations

The need for buffers in biological and biochemical research is universal. However, in the past, very few buffers in the important pH range of 6 to 8 were available. Those that were available were inappropriate for biological research and had serious disadvantages, such as toxicity or undesired reactivity. Phosphate buffers, for example, exhibit poor buffering capacity above pH 7.5, and they often inhibit reactions and precipitate polyvalent cations. Below pH 7.5, buffers such as TRIS can be toxic and show poor buffering capacity. Similarly, glycylglycine is useful above pH 8, but is of no value below pH 7.5.

Several important criteria must be met in order for a buffer to be useful in biological systems:

- 1 The buffers must be enzymatically and hydrolytically stable.
- 2 The pK_a of the buffer should be between 6 and 8 for most biological reactions.
- 3 The pH of the buffer solution should be minimally affected by concentration, temperature, ionic composition, or salt effects of the medium.
- 4 The buffer should be soluble in water and relatively insoluble in other solvents.
- 5 Cationic complexes should be soluble.
- 6 The buffer should exhibit no absorption of light in either the visible or UV regions.

Some years ago, Good¹ described a series of zwitterionic buffers possessing these characteristics. These so-called "Good's Buffers" are now widely used in cell culture, electrophoresis, biological systems and biochemical reactions. Over the years, several new zwitterionic buffers have been added to the original list of Good's buffers, and a list of these is shown in Table 1.

pK_a	Buffer	Cat. No.	pH Range	MW	Water Solubility (0°C, gm/100 mL)
6.15	MES	ICN195309	5.8 - 6.5	195.2	12.7
6.50	BIS-TRIS	ICN101038	5.8 - 7.2	209.2	20.9
6.76	PIPES	ICN190257	6.1 - 7.5	302.4	slightly
7.15	BES	ICN10092705, 5 g	6.6 - 7.6	213.2	68.2
		ICN10092780, 100 g			
		ICN10092783, 250 g			
		ICN10092791, 1 kg			
		MP210092725, 25 g			
7.20	MOPS	ICN102370	6.5 - 7.9	209.3	6.5
7.55	HEPES	ICN101926	7.0 - 8.0	238.3	53.6
7.80	HEPPSO	ICN151236	7.1 - 8.5	268.3	26.8
8.00	HEPPS	ICN10192725, 25 g	7.6 - 8.6	252.3	39.9
8.10	TRIS	MP21521761, 100 g	7.0 - 9.0	121.1	50.0
		MP21521765, 500g			
		ICN15217601, 1 kg			
		ICN15217605, 5 kg			
		ICN15217610, 10 kg			
8.15	TRICINE	ICN103112	7.6 - 8.8	179.2	14.3
8.35	BICINE	ICN10100525, 25 g	7.8 - 8.8	163.2	18.0
		ICN10100580, 100 g			
		ICN10100583, 250 g			
		MP210100591, 1 kg			
10.40	CAPS	ICN101435	9.7 - 11.1	221.3	10.4

Table 1. Biological and Biochemical Buffers

Zwitterionic buffers are typically supplied in the free acid form, although several are available as sodium salts, to aid in their solubility. As a general rule, a buffer is chosen so that the pK_a is slightly below the desired pH. By then adjusting with a suitable base, the buffer is brought to the desired pH.

Tissue Culture Applications

Several of the Good's buffers, most notably HEPES, TRICINE and TES, have been shown to be very effective in cell culture. Ceccarini and Eagle² have studied the optimum pH for growth of a number of normal, virus-transformed, and cancer cells, using various zwitterionic buffers to stabilize pH.

A study by Eagle³ has shown that eight of the Good's buffers are non-toxic. These buffers include BIS-TRIS, PIPES, BES, TES, HEPES, HEPPS, TRICINE and Bicine. A table of suggested buffer combinations for use in the presence of bicarbonate is also presented in Eagle's study.

In a study by Shipman⁴, HEPES was found to give higher maximum cell densities and viabilities in cultures, such as human embryonic lung, chick embryo fibroblast and guinea pig spleen cells. In viral studies, Shipman also observed that HEPES-buffered saline did not affect Rubella virus titration or hemagglutination assays for Polyoma or Sendai viruses. Phosphate-buffered saline had been reported to affect these determinations.

Description	CAS #	Formula	MW	Size	Cat. No.
BES [N,N-bis(2-Hydroxyethyl)-2-aminoethanesulfonic acid]. Free Acid. $pK_a = 7.15$. Useful pH range 6.6–7.6. BES buffer has been used in calcium phosphate-mediated transfection of eukaryotic cells with plasmid DNA.	[10191-18-1]	$C_6H_{15}NO_5S$	213.3	5 g	ICN10092705
				25 g	MP210092725
				100 g	ICN10092780
				250 g	ICN10092783
				1 kg	ICN10092791
BICINE [N,N-bis(2-Hydroxyethyl)glycine]. $pK_a = 8.35$. Useful pH range 7.8–8.8. BICINE is used in protein crystallization, studying enzyme reactions and electrophoresis.	[150-25-4]	$C_6H_{13}NO_4$	163.2	25 g	ICN10100525
				100 g	ICN10100580
				500 g	ICN10100583
				1 kg	MP210100591
BIS-TRIS [2,2-bis(Hydroxymethyl)-2,2',2''-nitrilotriethanol]. $pK_a = 6.50$. Useful pH range 5.8–7.2. A zwitterionic buffer used to calibrate glass electrodes and for nucleic acid and protein crystallizations.	[6976-37-0]	$C_8H_{19}NO_5$	209.2	25 g	ICN10103825
				100 g	ICN10103880
				500 g	ICN10103890
				1 kg	ICN10103891
CAPS [3-(Cyclohexylamino)propanesulfonic acid]. $pK_a = 10.4$. Useful pH range 9.7–11.1. A zwitterionic buffer used for protein sequencing and identification, Western blotting and immunoblotting.	[1135-40-6]	$C_9H_{19}NO_3S$	221.3	25 g	ICN10143525
				100 g	ICN10143580
				1 kg	ICN10143591
HEPES (N-2-Hydroxyethyl)piperazine-N'-2-ethanesulfonic acid). $pK_a = 7.55$. Useful pH range 7.0–8.0. A zwitterionic Good's buffer widely used in cell culture media and as an ampholytic separator to create a pH gradient in isoelectric focusing.	[7365-45-9]	$C_8H_{18}N_2O_4S$	238.3	25 g	ICN10192625
				100 g	ICN10192680
				250 g	ICN10192683
				1 kg	ICN10192691

Practical Considerations

Description	CAS #	Formula	MW	Size	Cat. No.
<p>HEPES HEMISODIUM SALT (N-2-Hydroxyethylpiperazine-N'-2-ethanesulfonic acid). Hemisodium salt. $pK_a = 7.5$. Useful pH range 6.8–7.2. Zwitterionic buffer widely used to maintain physiological pH, with slightly better solubility than HEPES free acid.</p>	[103404-87-1]	$C_8H_{17}N_2O_4 \cdot 1/2Na$	249.3	100 g	ICN15245180
<p>HEPES SODIUM SALT (N-2-Hydroxyethylpiperazine-N'-3-ethanesulfonic acid). Sodium salt. $pK_a = 7.5$. Useful pH range 6.8–7.2. Zwitterionic buffer widely used to maintain physiological pH, with slightly better solubility than HEPES free acid.</p>	[75277-39-3]	$C_8H_{17}N_2O_4Na$	260.3	25 g 100 g 250 g 1 kg	ICN10559325 ICN10559380 ICN10559383 ICN10559391
<p>HEPPS (N-2-Hydroxyethylpiperazine-N'-3-propanesulfonic acid). $pK_a = 8.00$. Useful pH range 7.6–8.6. This is the propane analog of HEPES and has many similar properties. Suitable for use in phosphorylation reactions when metal binding may occur. In mice it has been shown to break-up amyloid beta plaques associated with Alzheimer's Disease.</p>	[16052-06-5]	$C_9H_{20}N_2O_4S$	252.3	25 g	ICN10192725
<p>HEPPSO [4-(2-Hydroxyethyl)piperazine-1-(2-hydroxypropanesulfonic acid)]. $pK_a = 7.80$. Useful pH range 7.1–8.5. Zwitterionic buffer commonly used as an ampholytic separator to create a pH gradient in isoelectric focusing.</p>	[68399-78-0]	$C_9H_{20}N_2O_5S$	268.3	10 g 25 g 100 g	ICN15123610 ICN15123650 ICN15123680
<p>MES [2-(N-Morpholino)ethanesulfonic acid]. Monohydrate. $pK_a = 6.15$. Useful pH range 5.8–6.5. A zwitterionic buffer used in SDS-PAGE applications, preparation of culture media, and fluorescence microscopy. One of the first Good's buffers used for protein purification.</p>	[4432-31-9]	$C_6H_{13}NO_4S \cdot H_2O$	213.2	25 g 100 g 250 g 1 kg	ICN19530925 ICN19530980 ICN19530983 ICN19530991
<p>MES SODIUM SALT [2-(N-Morpholino)ethanesulfonic acid]. Sodium salt. $pK_a = 6.15$. Useful pH range 5.8–6.5. A zwitterionic buffer used in SDS-PAGE applications, preparation of culture media, and fluorescence microscopy. One of the first Good's buffers used for protein purification.</p>	[71119-23-8]	$C_6H_{12}NO_4SNa$	217.2	10 g 100 g	ICN15245410 ICN15245480
<p>MOPS [3-(N-Morpholino)propanesulfonic acid]. Free Acid. $pK_a = 7.20$. Useful pH range 6.5–7.9. Widely used zwitterionic buffer due to its inert properties. It does not interact with any metal ions in solution. Used in mammalian cell culture and denaturing gel electrophoresis of RNA. Interacts with BSA and stabilizes it.</p>	[1132-61-2]	$C_7H_{15}NO_4S$	209.3	100 g 250 g 1 kg	ICN10237080 ICN10237083 ICN10237091

Description	CAS #	Formula	MW	Size	Cat. No.
<p>MOPS SODIUM SALT [3-(N-Morpholino)propanesulfonic acid]. Sodium Salt. pK_a = 7.20. Useful pH range 6.5–7.9. Widely used zwitterionic buffer in cell culture. MOPS can modify lipid interactions and influence the thickness and barrier properties of membranes. Interacts with BSA and stabilizes it.</p>	[71119-22-7]	C ₇ H ₁₄ NO ₄ SNa	231.2	25 g 100 g 250 g 1 kg	ICN19067025 ICN19067080 ICN19067083 ICN19067091
<p>MOPSO SODIUM SALT [3-(N-Morpholino)-2-hydroxypropane sulfonic acid]. Sodium Salt. pK_a = 6.90. Useful pH range 6.2–7.6. A zwitterionic buffer commonly used for cell culture media, as a running buffer in electrophoresis, and for protein purification. Although MOPSO does not form complexes with most metals, it may have a strong interaction with iron in solution.</p>	[79803-73-9]	C ₇ H ₁₄ NO ₅ SNa	247.2	100 g	ICN15245580
<p>PIPES [Piperazine-N,N'-bis(2-ethanesulfonic acid)]. Free Acid. pK_a = 6.76. Useful pH range 6.1–7.5. A zwitterionic buffer used in cell culture and protein purification. PIPES can minimize lipid loss when buffering glutaraldehyde histology in plant and animal tissues.</p>	[5625-37-6]	C ₈ H ₁₈ N ₂ O ₆ S ₂	302.4	100 g 500 g 1 kg	ICN19025780 ICN19025790 ICN19025791
<p>TRICINE [N-tris(Hydroxymethyl)methylglycine]. pK_a = 8.15. Useful pH range 7.6–8.8. A zwitterionic buffer used in SDS-PAGE procedures to separate low molecular weight peptides.</p>	[5704-04-1]	C ₆ H ₁₃ NO ₅	179.2	25 g 100 g 250 g 1 kg	ICN10311225 ICN10311280 ICN10311283 ICN10311291
<p>TRIS [Tris-(hydroxymethyl)aminomethane; Tromethamine; Trometamol]. Purity: 99.0–99.5%. pK_a = 8.1. Useful pH range 7.0–9.0. Widely used buffer component for buffer solutions and protein purification. This grade of TRIS is excellent where purity and value are both important. It is superior to technical grade and less expensive than Ultra Pure material.</p>	[77-86-1]	C ₄ H ₁₁ NO ₃	121.1	100 g 500 g 1 kg 5 kg	MP21521761 MP21521765 ICN15217601 ICN15217605
<p>TRIS USP [Tris-(hydroxymethyl)aminomethane; Tromethamine; Trometamol]. USP Grade. Purity: 99.95% minimum. pK_a = 8.1. Useful pH range 7.0–9.0. Excellent biochemical and biological buffer where certified high purity is required.</p>	[77-86-1]	C ₄ H ₁₁ NO ₃	121.1	100 g 500 g 1 kg	ICN19560580 ICN19560590 ICN19560591
<p>TRIS ULTRA PURE [Tris-(hydroxymethyl)aminomethane; Tromethamine; Trometamol]. Ultra Pure Grade. Purity: 99.95% minimum. pK_a = 8.1. Useful pH range 7.0–9.0. Excellent biochemical and biological buffer for all applications where high purity is required.</p>	[77-86-1]	C ₄ H ₁₁ NO ₃	121.1	100 g 500 g 1 kg	MP21031331 MP21031335 ICN10313301

References

- Roy, R. N.; et al. Buffer Standards for the Physiological pH of the Zwitterionic Compound, ACES, from 5 to 55 °C. *Journal of Solution Chemistry*. **2009**, *38*, 471-483.
- Zawisza, I.; et al. Cu(II) complex formation by ACES buffer. *Journal of Inorganic Biochemistry*. **2013**, *129*, 58-61.
- Huhta, E.; Parjanen, A.; Mikkola, S. A kinetic study on the chemical cleavage of nucleoside diphosphate sugars. *Carbohydrate Research*. **2010**, *345*, 696-703.
- Williams, T. I.; et al. A novel BICINE running buffer system for doubled sodium dodecyl sulfate – polyacrylamide gel electrophoresis of membrane proteins. *Electrophoresis*. **2006**, *27* (14).
- Newman, J. Novel buffer systems for macromolecular crystallization. *Acta Cryst*. **2004**, *D60*, 610-612.
- Goklen, K. E.; Suda, E. J.; Ubiera, A. R. Buffer system for protein purification. US Patent 9,624,261, 2017.
- Naz, H.; et al. Effect of pH on the structure, function, and stability of human calcium/calmodulin-dependent protein kinase IV: combined spectroscopic and MD simulation studies. *Biochemistry and Cell Biology*. **2016**, *94*(3), 221-228.
- Lee, S. Y.; et al. Densities, Viscosities, and Refractive Indexes of Good's Buffer Ionic Liquids. *J. Chem. Eng. Data*. **2016**, *61*, 2260-2268.
- Lopez, A.; Liu, J. DNA-templated fluorescent gold nanoclusters reduced by Good's buffer: from blue-emitting seeds to red and near infrared emitters. *Canadian Journal of Chemistry*. **2015**, *93*(6), 615-620.
- Wang, D.; et al. Stability study of tubular DNA origami in the presence of protein crystallisation buffer. *RSC Adv*. **2015**, *5*, 58734-58737.
- Schmidt, J.; et al. Effect of Tris, MOPS, and phosphate buffers on the hydrolysis of polyethylene terephthalate films by polyester hydrolases. *FEBS Open Bio*. **2016**, *6* (9), 919-927.
- Dias, K. C.; et al. Influence of different buffers (HEPES/MOPS) on keratinocyte cell viability and microbial growth. *Journal of Microbiological Methods*. **2016**, *125*, 40-42.
- Chandra, K.; et al. Separation of Stabilized MOPS Gold Nanostars by Density Gradient Centrifugation. *ACS Omega*. **2017**, *2* (8), 4878-4884.
- Grabber, J. H. Relationships between Cell Wall Digestibility and Lignin Content as Influenced by Lignin Type and Analysis Method. *Crop Sci*. **2019**, *59*, 1122-1132.
- Ahmed, S. R.; et al. Synthesis of Gold Nanoparticles with Buffer-Dependent Variations of Size and Morphology in Biological Buffers. *Nanoscale Research Letters*. **2016**, *11*, 65.
- Pannuru, P.; et al. The effects of biological buffers TRIS, TAPS, TES on the stability of lysozyme. *International Journal of Biological Macromolecules*. **2018**, *12*, 720-727.
- Samuelsen, L.; et al. Buffer solutions in drug formulation and processing: How pKa values depend on temperature, pressure and ionic strength. *International Journal of Pharmaceutics*. **2019**, *560*, 357-364.
- Venturini, E.; et al. Targeting the Potassium Channel Kv1.3 Kills Glioblastoma Cells. *Neurosignals*. **2017**, *25*, 26-38.
- Haider, S. R.; Reid, H. J.; Sharp, B. L. Tricine-SDS-PAGE. In *Electrophoretic Separation of Proteins*; Kurien, B., Scofield, R., Eds.; Methods in Molecular Biology; Humana Press: New York, 2019; Vol. 1855.
- Ying, Y.; et al. Solubilization of proteins in extracted oil bodies by SDS: A simple and efficient protein sample preparation method for Tricine-SDS-PAGE. *Food Chemistry*. **2015**, *181*, 179-185.
- Taha, M.; et al. Good's buffer ionic liquids as relevant phase-forming components of self-buffered aqueous biphasic systems. *Journal of Chemical Technology & Biotechnology*. **2017**, *92* (9), 2287-2299.
- Gao, Y.; et al. Highly effective electrochemical water oxidation by copper oxide film generated in situ from Cu(II) tricine complex. *Chinese Journal of Catalysis*. **2018**, *39* (3), 479-486.
- Taha, M.; Khan, I.; Coutinho, J. A. P. Coordination abilities of Good's buffer ionic liquids toward europium(III) ion in aqueous solution. *Journal of Chemical Thermodynamics*. **2016**, *94*, 152-159.
- Takeshita, Y.; et al. The effects of pressure on pH of Tris buffer in synthetic seawater. *Marine Chemistry*. **2017**, *188*, 1-5.
- Ibrahim-Hashim, A.; et al. Tris-base buffer: a promising new inhibitor for cancer progression and metastasis. *Cancer Medicine*. **2017**, *6* (7), 1720-1729.
- Anes, B.; Bettencourt da Silva, R. J. N.; Oliveira, C.; Camões, M. F. Seawater pH measurements with a combination glass electrode and high ionic strength TRIS-TRIS HCl reference buffers – An uncertainty evaluation approach. *Talanta*. **2019**, *193*, 118-122.

Ultra Pure Reagents

For critical, sensitive, demanding research where even a very minute amount of contaminant can potentially wreak havoc, MP Biomedicals Ultra Pure Reagents can provide the high quality you require. Using special purification steps, such as multiple re-distillations and recrystallizations (up to 5X), MP Bio purifies these reagents to uncommonly stringent specifications, making these products truly Ultra Pure. For example, during gel electrophoresis, it is often difficult to work at lower temperatures and pH because of marked precipitation when using sodium dodecyl sulfate (SDS). MP Bio solves this problem with our Ultra Pure lithium dodecyl sulfate (LDS), which exhibits greater solubility than SDS at lower temperatures, while maintaining similar detergency and wetting ability. Substitution of Ultra Pure LDS for SDS has been shown to result in greater resolution for certain proteins. Similarly, metallic and anionic contaminants, even in minute amounts, can shut down or block enzymatic proteins, resulting in poor yields and/or incorrect analytical and electrophoretic results. Use of Ultra Pure reagents often eliminates trace amounts of metallic contaminants and provides a better result. Remember, if it doesn't say "Ultra Pure", it probably isn't. With MP Bio Ultra Pure reagents, no finer quality products are available anywhere, at any price.

Name	Description	Pack Size	Cat. No.
Acrylamide, Ultra Pure	C ₃ H ₅ NO MW 71.1. Purity >99.9%. Acrylic acid content: < 0.001%. Super pure monomer for preparation of polyacrylamide gels for sensitive PAGE applications.	500 g	MP04814326
Ammonium Sulfate, Ultra Pure	(NH ₄) ₂ SO ₄ MW 132.2. Purity: ≥ 99%. A widely used reagent in molecular biology for the isolation and purification of enzymes and proteins. It is used for the precipitation or fractionation of proteins and for purification of antibodies. Ammonium sulfate is used in long PCR buffer, in PCR lysis solution, and in second strand cDNA synthesis.	50 g	MP04808210
		1 kg	ICN808211
		5 kg	MP04808229
Cesium Chloride, Ultra Pure	CsCl MW 168.36. Purity: ≥99.999%. Cesium chloride is typically used for density gradient work and for the purification of virus/phage, nucleic acids and nucleoproteins. It is used for the preparation of electrically conducting glasses, used to make solutions for the separation of RNA from DNA by density gradient centrifugation.	5 g	MP215058905
		25 g	MP215058925
		100 g	MP215058980
		500 g	MP215058990
		1 kg	MP215058991
Guanidine Hydrochloride, Ultra Pure	Purity: ≥ 99.5%. This strong denaturant can solubilize insoluble or denatured proteins, such as inclusion bodies. Highly concentrated (6 - 8 M) Guanidine HCl solutions are used to denature native globular proteins, presumably by disrupting the hydrogen bonds that hold the protein in its unique structure.	25 g	ICN10569625
		100 g	ICN10569680
		500 g	ICN10569690
		1 kg	ICN10569691
		5 kg	ICN10569694
N-Lauroylsarcosine sodium salt, Ultra Pure	Purity: ≥97%. An anionic detergent useful in the cell lysis process of RNA purification. Ideal for solubilizing membrane proteins prior to electrophoresis.	50 g	ICN19400950
		100 g	ICN19400980
		500 g	ICN19400990
Lithium dodecylsulfate, Ultra Pure	(LDS). Purity: >99%. Detergent for solubilizing proteins for electrophoresis. Demonstrates greater solubility than SDS at lower temperatures, while maintaining similar detergency and wetting ability.	5 g	ICN800752
		25 g	ICN800753

Ultra Pure Reagents

Name	Description	Pack Size	Cat. No.
N,N'-Methylene-bis-acrylamide, Ultra Pure	Purity: 99.9%. A highly purified bisacrylamide for crosslinking with acrylamide to make superior PAGE gels for critical electrophoresis applications. May be used in UV scanning gels due to its optical clarity. Acrylic acid content: <0.02%	5 g	ICN800172
Phenol, Ultra Pure, 99%	For the extraction of nucleic acids and to solubilize and denature proteins. Typically used in a mixture of phenol and buffered aqueous solution, proteins are denatured and collected at the interphase, while most nucleic acids remain in the aqueous phase.	500 g	MP04818048
		1 kg	ICN800673
Sodium dodecylsulfate, Ultra Pure	Purity: ≥99%. An anionic surfactant that denatures and solubilizes proteins for electrophoresis. Also useful as an aid in cell lysis during DNA extraction, and for dispersing and suspending nanotubes.	25 g	MP04811033
		50 g	ICN811036
		100 g	ICN811034
		500 g	MP04811032
		1 kg	MP04811030
Sucrose, Ultra Pure	C ₁₂ H ₂₂ O ₁₁ M.W. 342.30. Purity: 99.9%. DNase and RNase-free. Used for preparation of sucrose gradients for purification of proteins and RNAs.	100 g	NC1637589
		500 g	ICN821713
		1 kg	MP04821721
Tris(hydroxymethyl) aminomethane, Ultra Pure, 99.95%	(TRIS base). Purity: 99.95%. Widely used zwitterionic Good's buffer for preparation of many different electrophoresis buffers. pK _a = 8.06 at 20°C.	100 g	MP21031331
		250 g	MP21031332
		500 g	MP21031335
		1 kg	ICN10313301
		5 kg	ICN10313305
Tris(hydroxymethyl) aminomethane, Ultra Pure, 99.9%	(TRIS base). Purity: 99.9%. Widely used zwitterionic Good's buffer for preparation of many different electrophoresis buffers. pK _a = 8.06 at 20°C.	500 g	MP04819620
		1 kg	MP04819623
		5 kg	ICN819638
Urea, Ultrapure, 99%	Purity: 99%. A high purity protein denaturant frequently added to buffers and solutions used in protein research.	1 lb	ICN10569501
		5 lb	ICN10569505
Urea, Ultra Pure	CH ₄ N ₂ O M.W. 60.06. Purity: ≥99%. An ultra pure reagent suitable for use as a protein denaturant. Urea is commonly used to solubilize and denature proteins for denaturing isoelectric focusing and two-dimensional electrophoresis and in acetic acid-urea PAGE gels. Urea is typically used at a concentration of 8 M for protein denaturation or solubilization. A final concentration of 5 M urea is commonly used in molecular biology for sequencing gels.	1 lb	MP04821519
		5 lb	ICN821527
		25 lb	MP04821532
		1 kg	MP04821528
		5 kg	ICN821530
		10 kg	ICN821858

Make Your Own Buffer Solutions Cookbook

The following are recommended recipes for preparing the most commonly used buffers in electrophoresis applications. Whenever possible, MP Bio strongly recommends using Ultra Pure reagents and water when preparing them.

Tris-Glycine Native Running Buffer

Format: **500 mL of 10X solution** Shelf-life: **1 year at room temperature** pH: **8.3**

Component	1X Concentration	Quantity for 10X solution
Tris	25 mM	29.0 g
Glycine	192 mM	144.0 g
Deionized water (ultra pure)	—	to 1.0 L

Tris-Glycine Native Sample Buffer

Format: **20 mL of 2X solution** Shelf-life: **1 year at 4°C** pH: **8.6**

Component	1X Concentration	Quantity for 2X solution
Tris HCl	100 mM	4 mL of a 0.5 M sol.
Glycerol	10%	2 mL
Bromophenol Blue	0.0025%	0.5 mL of a 1% sol
Deionized water (ultra pure)	—	to 10.0 mL

Tris-Glycine Native Transfer Buffer

Format: **500 mL of 25X solution** Shelf-life: **1 year at room temperature** pH: **8.3**

Component	1X Concentration	Quantity for 25X solution
Tris	12 mM	18.2 g
Glycine	96 mM	90.0 g
Deionized water (ultra pure)	—	to 500 mL

Tris-Glycine-SDS Running Buffer

Format: **500 mL of 10X solution** Shelf-life: **1 year at room temperature** pH: **8.3**

Component	1X Concentration	Quantity for 10X solution
Tris	25 mM	29.0 g
Glycine	192 mM	144.0 g
SDS	0.1%	10.0 g
Deionized water (ultra pure)	—	to 1.0 L

Tris-Glycine-SDS Sample Buffer

Format: **20 mL of 2X solution** Shelf-life: **1 year at 4°C** pH: **6.8**

Component	1X Concentration	Quantity for 2X solution
Tris HCl	63 mM	2.5 mL of a 0.5 M sol.
Glycerol	10%	2 mL
SDS	2%	4 mL of a 10% (wv) Sol.
Bromophenol Blue	0.0025%	0.5 mL of a 1% Sol.
Deionized water (ultra pure)	—	to 10.0 mL

Tris-Tricine-SDS Running Buffer

Format: **500 mL of 10X solution** Shelf-life: **1 year at room temperature** pH: **8.3**

Component	1X Concentration	Quantity for 10X solution
Tris pH 8.3	100 mM	121.0 g
Tricine	100 mM	179.0 g
SDS	0.1%	10.0 g
Deionized water (ultra pure)	—	to 1.0 L

Tris-Tricine-SDS Sample Buffer

Format: **20 mL of 2X solution** Shelf-life: **1 year at 4°C** pH: **8.45**

Component	1X Concentration	Quantity for 2X solution
Tris HCl, pH 8.45	450 mM	3 mL of a 3.0 M sol.
Glycerol	12%	2.4 mL
SDS	4%	0.8 g
Coomassie Blue G250	0.0025%	0.5 mL of a 1% sol.
Phenol Red	0.0025%	0.5 mL of a 1% sol.
Deionized water (pure water)	—	to 10.0 mL

TBE Running Buffer

Format: **1000 mL of 5X solution** Shelf-life: **1 year at room temperature** pH: **8.3**

Component	1X Concentration	Quantity for 5X solution
Tris	89 mM	54.0 g
Boric acid	89 mM	27.5 g
EDTA (free acid)	2 mM	2.9 g
Deionized water (ultra pure)	—	to 1.0 L

TBE Sample Buffer

Format: **10 mL of 6X solution** Shelf-life: **1 year at 4°C**

Component	1X Concentration	Quantity for 6X solution
Tris	45 mM	6 mL of 5X TBE running buffer
Boric acid	45 mM	—
EDTA (free acid)	1 mM	—
Glycerol	5.3%	3.2 mL
Bromophenol Blue	0.005%	0.3 mL of a 1% Sol.
Xylene Cyanol	0.005%	0.3 mL of a 1% Sol.
Deionized water (ultra pure)	—	to 10.0 mL



One Call. One Source. A World of Ultra Pure Biochemicals.

- Apoptosis
- Cell Biology
- Culture Growth Media
- FastPrep® Sample Prep
- Immunology
- Molecular Biology
- Adsorbents
- Biochemicals
- Fine Chemicals
- Labware
- Dosimetry
- Research Diets
- SafTest™ Food Quality
- Diagnostics
- Drugs of Abuse
- Infectious Disease
- EIA/RIA

Contact us today:

In the United States

Order online: fishersci.com
Fax an order: 1-800-926-1166
Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca
Fax an order: 1-800-463-2996
Call customer service: 1-800-234-7437

