



## New Millicell<sup>®</sup> EZ SLIDE

- Culture directly on a microscope slide
- Simplify histology and cell staining procedures
- Simple detachment of media chamber from slide



**Simplify your cell analysis** by using the Millicell EZ SLIDE to culture, fix, stain and view your sample all in one device. There's no need to remove the media chamber from the slide prior to fixing or staining. Easy well removal eliminates the worry of breaking slides or harming cells. Acquire data simply and quickly with Millicell EZ SLIDES.

### Features

- Breakable tabs for easy well removal
- No tools required to remove wells from slide
- No adhesive used to attach wells to slide
- No leaky wells
- Polypropylene wells allow for fixing and staining with wells still attached
- Slide holder protects slide from damage

*This new, disposable, molded, culture slide assembly is compatible with all growth media, and most fixing solutions and staining reagents.*


### Free Samples Available!

Contact your local Millipore sales representative or visit [www.millipore.com/EZSLIDE](http://www.millipore.com/EZSLIDE)

(See reverse for contact information.)

## Easily Remove Wells Before Microscopy:

- 1 Hold the Millicell EZ SLIDE down on a flat surface.
- 2 Break away four tabs by pushing them down with your thumb.
- 3 Remove the cover.
- 4 Hold the slide down and lift the wells upward to remove.
- 5 Remove slide from holder.



## ORDERING INFORMATION

Description	Qty/Pk	Catalogue No.
Millicell EZ SLIDE (4-well glass)	16	PEZGS0416
	96	PEZGS0496
Millicell EZ SLIDE (8-well glass)	16	PEZGS0816
	96	PEZGS0896
<b>Accessory</b>		
Millicell EZ SLIDE Microscope Slide Holder	1	PEZXMSH01

[www.millipore.com/EZSLIDE](http://www.millipore.com/EZSLIDE)



For technical assistance, contact Millipore:  
**1-800-MILLIPORE (1-800-645-5476)**  
E-mail: [tech\\_service@millipore.com](mailto:tech_service@millipore.com)



For customer service, call **1-800-766-7000**.  
To fax an order, use **1-800-926-1166**.  
To order online: [www.fishersci.com](http://www.fishersci.com)