




Catch the BUZZ

NEW! TransIT®-Insect Transfection Reagent

Insect cell expression is a convenient vehicle for high volume mammalian protein production. Recombinant baculovirus production and transient transfection are commonly used methods for expression in insect cells.

TransIT-Insect Transfection Reagent offers distinct advantages over other reagents such as:

- Exceptional DNA Delivery In insect cell types including Sf9, High FiveTM and S2
- High Baculovirus Production Ideal for baculovirus expression in insect cells
- Serum Compatibility Non-liposomal, animal-origin free formulation
- Better Value Low reagent amounts required per transfection




Efficient Delivery of Baculovirus Genomic DNA. Sf9 cells were co-transfected with $0.5~\mu g$ of ProGreenTM baculovirus genomic vector DNA (AB Vector) encoding green-fluorescent protein (GFP) and $0.1~\mu g$ of pVL1393 transfer vector (AB Vector) using $TransT^{\odot}$ -Insect Transfection Reagent at the reagent-to-total DNA ratio of $3.1~(\mu l: \mu g)$.



ADVANCE YOUR TRANSFECTIONS.


Request a FREE SAMPLE of *Trans*IT®-Insect Transfection Reagent. Contact your local Fisher Representative


Mirus Bio - The Transfection Experts.





Superior Gene Expression Compared to Liposomal Formulations





TransIT®-Insect Outperforms Competitor Transfection Reagents. Insect cell lines Sf9, High Five™, and *Drosophila melanogaster* Schneiders 2 (S2 or D. Mel. (2)) were transfected with a luciferase expression plasmid driven by an hr5 enhancer/IE1 promoter using the designated reagent at the indicated reagent-to-DNA ratios. Transfections were performed in 96-well plates using 0.1 µg of plasmid DNA. Luciferase expression was measured at 48 hours post-transfection using a standard assay. Sf9 and High Five cells were cultured and transfected in serum-free media formulations; S2 cells were cultured and transfected in serum containing medium. Error bars represent the standard error of the mean for triplicate wells.

Insect Cell DNA Transfection

Product Name	Product No.	Quantity
TransIT®-Insect	MIR6104	0.4 ml
Transfection Reagent	MIR6100	1 ml
	MIR6105	5 x 1 ml
	MIR6106	10 x 1 ml

To inquire about bulk pricing, please call 888.530.0801

FLTIInsect-12/2013 BN1209132

©2013. All rights reserved. Mirus Bio LLC. All trademarks are the property of their respective owners. TransIT is a registered trademark of Mirus Bio LLC. High Five is a trademark and Cellfectin is a registered trademark of Life Technologies Corporation. ProGreen is a trademark of AB Vector LLC. ESCORT is a trademark of Sigma Aldrich Co LLC.

Prove it to Yourself with a FREE Sample



Start with our Reagent Agent® to determine the best solution for your transfection needs: mirusbio.com/ra

Mirus Bio products are also distributed in the U.S. by:



For customer service, call 1-800-766-7000. To fax an order, use 1-800-926-1166. To order online: www.fishersci.com