BD Biosciences Fluorochrome Reference Chart

Fluorescence

Excitation Laser

Visit bdbiosciences.com/colors for detailed information about our newest fluorochromes and instrumentation.

To select your optimal combination of fluorochromes, visit bdbiosciences.com/spectra to use an interactive fluorescence spectrum tool.

Instrument	Laser	Line (nm) Fluorescence Channel		Fluorochromes provided by BD Biosciences		
BD FACSArray™	Green Diode	532	Yellow	PE		•
bioanalyzer	Red Diode	635	Far Red Red	PerCP-Cy5.5 APC	PE-Cy7 Alexa Fluor® 647	
	Ned Diode	055	Infrared	BD APC-H7	APC-Cy7	
*BD FACSCalibur™	Argon	488	FL1 Green	FITC	Alexa Fluor® 488	
flow cytometry system			FL2 Yellow FL3 Red	PE-Cy5 ^a	PerCP	PerCP-Cy5.5 PE-Cy7
	Red Diode	635	FL4 Red	APC ^a	Alexa Fluor® 647	Teleficitys.5
*BD FACSCanto™	Solid State	488	Green	FITC	Alexa Fluor® 488	
flow cytometry system			Yellow Red	PE PerCP	PerCP-Cy5.5	
			Infrared	PE-Cy7	rerer eysis	
	HeNe	633	Red	APC	Alexa Fluor® 647	
*†BD FACSCanto™ II	Solid State	488	Infrared Green	BD APC-H7 FITC	APC-Cy7 Alexa Fluor® 488	
flow cytometry system	John State	100	Yellow	PE	Alexa Hadre 100	
			Orange Red	PE-Texas Red®b	DowCD Cut t	•
			Infrared	PerCP PE-Cy7	PerCP-Cy5.5	
	HeNe	633	Red	APC	Alexa Fluor® 647	
			Far Red Infrared	Alexa Fluor® 700 ^b BD APC-H7	APC-Cy7	
	Solid State ^b	405	Green	BD Horizon [™] 500 ^b	AmCyan ^b	
			Blue	BD Horizon™ V450b	Pacific Blue ^{™,b}	
Preconfigured BD TM LSR II	Solid State	488	Green Yellow	FITC	Alexa Fluor® 488	
(typical setup) ^d			Orange	PE-Texas Red®		
			Red	PerCP	PE-Cy5ª	PerCP-Cy5.5
	Solid State	640	Infrared Red	PE-Cy7 APC ^a	Alexa Fluor® 647	
	Solid State	040	Far Red	Alexa Fluor® 700	AICAG FIGOR 647	•
			Infrared	BD APC-H7	APC-Cy7	
	Solid State	405	Green Blue	BD Horizon V500 BD Horizon V450	AmCyan Pacific Blue™	
Special Order BD™ LSR II	Solid State	488	Green	FITC	Alexa Fluor® 488	
Special Order			Yellow	PE D. LO		
BD LSRFortessa TM			Orange Red	PE-Texas Red® PerCP	PE-Cy5ª	PerCP-Cy5.5
(typical setup) ^d		_	Infrared	PE-Cy7	Í	,
	Solid State	532 or 561	Yellow	PE PE-Texas Red®		
			Orange Red	PE-Cy5 ^a		
			Infrared	PE-Cy7		
	Solid State	640	Red Far Red	APC ^a Alexa Fluor® 700	Alexa Fluor® 647	
			Infrared	BD APC-H7	APC-Cy7	I
	Solid State	405	Green Blue	BD Horizon V500 BD Horizon V450	AmCyan Pacific Blue™	
BD FACSAria™ cell sorter	Solid State	488	Green	FITC	Alexa Fluor® 488	
family ^c			Yellow	PE		_
(typical setup) ^d			Orange Red	PE-Texas Red® PerCP	PE-Cy5 ^a	PerCP-Cy5.5
			Infrared	PE-Cy7	TE Cy3	reici cys.s
	Solid State ^b	561	Yellow	PE		
			Orange Red	PE-Texas Red® PE-Cy5ª		
			Infrared	PE-Cy7		
	HeNe	640	Red	APC ^a	Alexa Fluor® 647	
			Far Red Infrared	Alexa Fluor® 700 BD APC-H7	APC-Cy7	
	Solid State ^b	405	Green	BD Horizon V500	AmCyan	
BD Influx™ cell sorter	Solid State	488	Blue	BD Horizon V450	Pacific Blue™	
ייייייייייייייייייייייייייייייייייייי	John State	700	Green Yellow	PE PE	Alexa Fluor® 488	
			Orange	PE-Texas Red®		
			Red Infrared	PE-Cy5 PE-Cy7	PerCP-Cy5.5	
	Solid State	532 or 561	Yellow	PE		
			Orange	PE-Texas Red®		
			Red Infrared	PE-Cy5 PE-Cy7		
	Solid State	640	Red	APC	Alexa Fluor® 647	
			Far Red	Alexa Fluor®700	ADS 6.3	
	Solid State	405	Infrared Green	BD APC-H7 BD Horizon V500	APC-Cy7 AmCyan	
			Blue	BD Horizon V450	Pacific Blue™	
34.DC			tion have the Little of			TM J DD FACCA TM "

^aAPC and PE-Cy5 may be used together on instruments with cross-beam compensation. bAvailable through laser and/or detector options. bBD FACSAria™ and BD FACSAria™ ll ^dMore laser and detector options are available through the Special Order Research Products (SORP) program.

Choose a winning combination - Guidelines for selecting reagents for multicolor flow cytometry

1 The basics: Know your instrument Reagent selection starts with your instrument configuration. The lasers and detectors in your configuration dictate how well your cytometer can excite and measure a given fluorochrome, and whether you have enough detectors to read out a given combination of

7 Fluorochromes: Go for the bright Rank available dyes according to their intrinsic brightness on a particular instrument (when configured with a specified set of lasers and filters).

3 Minimize spillover As soon as cells are stained with multiple reagents, spectral overlap (or spillover) becomes an issue. The more colors you attempt to resolve on any particular cell, the more spillover impacts sensitivity. We use compensation, an adjustment applied to all colors, to correct for spillover. For example, a cell population fluorescing only in FITC will show no PE fluorescence, on average, but will likely exhibit more spread in the PE detector after

compensation than completely

unstained cells.

Colors and specificities: 4 Define winning combinations Once the fluorochromes to be used have been defined, you can begin to match antibody specificities to particular fluorochromes. Generally, reserve the brightest fluorochromes for dim antigens, and vice versa, but avoid spillover from bright cell populations into detectors

requiring high sensitivity for those

3 APC-Cy7, and to a lesser extent, PE-Cy7, can degrade in the presence of light, fixative, and elevated temperatures so that they emit in the parent dye detector (APC or PE). By minimizing the exposure of samples to light, heat, and formaldehydebased fixatives, this problem can be largely avoided. For more stable tandem dyes, BD now offers BD APC-H7 conjugated antibodies.

O Use controls (such as fluorescence-minus-one, or FMO) to validate your selected multicolor reagent cocktail. FMO controls help define the contribution of spillover to the background in a given detector, and are therefore useful in gauging the sensitivity of that detector in the context of a certain reagent cocktail.

Validation


For additional guidelines, visit **bdbiosciences.com/colors** to download the Application Note "Selecting Reagents for Multicolor Flow Cytometry."

Stain index of various fluorochrome conjugates on a BD™ LSR II

Reagent	Clone	Filter	Stain Index
PE	RPA-T4	575/26	305
APC ¹	RPA-T4	660/20	263
PE-Cy TM 5 ²	RPA-T4	695/40	198
Alexa Fluor® 647¹	RPA-T4	660/20	184
PE-Cy™7	RPA-T4	780/60	122
PerCP-Cy™5.5²	RPA-T4	695/40	99
Alexa Fluor® 488 ³	RPA-T4	530/30	68
BD Horizon™ V450 ⁵	RPA-T4	450/50	65
Alexa Fluor® 700	RPA-T4	720/40	64
Pacific Blue™,5	RPA-T4	450/50	63
FITC ³	RPA-T4	530/30	43
AmCyan ⁶	RPA-T4	525/50	37
APC-Cy7 ⁴	RPA-T4	780/60	36
PerCP ²	RPA-T4	695/40	30
BD Horizon™ V500 ⁶	RPA-T4	525/50	27
BD APC-H7 ⁴	RPA-T4	780/60	25

Freshly isolated lymphocytes, stained with anti-human CD4 antibodies conjugated with various fluorochromes run on a BD™ LSR II flow cytometer. This chart is meant as a guideline of relative stain indices of various fluorochromes. Observed relative stain indices may vary depending on instrument configurations and reagents used.

1, 2, 3, 4, 5, 6 Fluorochromes listed with the same superscript number are read in the same detector, and thus would not normally be used in combination.



Stain Index = D/W

Resolution sensitivity (the ability to resolve a dim positive signal from background) depends upon the difference between positive and background peak means (D) and the spread of the background peak (W). W, and W, represent background peaks with different spreads. The stain index is a metric that captures both of these factors.

* For In Vitro Diagnostic Use.

† Seven- and eight-color assays on this device are for Research Use Only.

Unless otherwise specified, all products are for Research Use Only.

Class I (1) laser product

APC-Cy7: US patent 5,714,386 Cy™ is a trademark of Amersham Biosciences Corp.

Pacific Blue is a trademark and Alexa Fluor and Texas Red are registered trademarks of Molecular Probes, Inc. Alexa Fluor®, Texas Red®, and Pacific Blue™ dye antibody conjugates are sold under license from Molecular Probes, Inc., for research use only, excluding use in combination with microarrays, or as analyte specific reagents. The Alexa Fluor® dyes (except for Alexa Fluor® 430) and Pacific Blue™ dye are covered by pending and issued patents.

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company.

23-11509-00

BN12170910

To order products through Fisher Scientific, please add BDB to the beginning of each part number listed.





For customer service, call 1-800-766-7000. To fax an order, use 1-800-926-1166. To order online: www.fishersci.com