



Features:

- High quality glassware at affordable prices
- Made from borosilicate 3.3 glass
- Superior chemical durability
- Autoclavable at 121 C
- Easy to read graduations
- Manufactured and calibrated to quality standards

Description	Cat. No.	Price
Beaker - 50 ml	S15437	2.60
Beaker - 100 ml	S15438	3.10
Beaker - 250 ml	S15439	3.10
Beaker - 500 ml	S15441	5.15
Beaker - 1000 ml	S15442	9.40

Complete Glassware Line



A full line of ELABGlass from Eisco Scientific is available through Fisher Science Education.