



# Fisher Scientific

**eISCO**  
Trust | Deliver | Learn

Premium Industrial Clamps

## Features:

- Zinc die cast construction made from Zamac 3 Alloy + Aluminum Alloy
- Connecting rods stainless steel
- Complete non ferrous construction
- High grade non toxic powdered coating

## Benefit

- Won't rust and will last for years
- Can be used in tough scientific environments
- High strength, supports heavy loads

Part #	Description
04999132	Premium Three Fingers Clamp - Vinyl Coated
04999133	Premium Three Finger Clamp w/ Boss head - Vinyl Coated
04999140	Premium, Large Three Fingers Clamp - Vinyl Coated
04999141	Premium Large Three Finger Clamp w/ Boss head - Vinyl Coated
04999144	Premium Three Finger Clamp - Cork Lined
04999145	Premium Three Finger Clamp w/ boss head - Cork Lined
04999138	Premium Three Finger Clamp w/ dual screw adjustments - Vinyl Coated
04999139	Premium Three Finger Clamp w/ dual screw adjustments w/ boss head - Vinyl Coated
04999157	Premium Four Prong Clamp - Cork Lined
04999134	Premium Four Prong Clamp w/ boss head - Cork Lined
04999158	Premium Four Prong Clamp - Vinyl Coated
04999135	Premium Four Prong Clamp - Vinyl Coated w/ boss head
04999159	Premium Two Prong Clamp - Vinyl Coated w/ dual screw adjustments
04999136	Premium Two Prong Clamp w/ dual screw adjust and swivel boss head - Vinyl Coated
04999160	Premium Two Prong Clamp - Vinyl Coated
04999137	Premium Two Prong Clamp w/ boss head attachment - Vinyl Coated
04999142	Premium Two Prong Clamp - Cork lined
04999143	Premium Two Prong Clamp cork lined w/ boss head attachment - Cork lined
04999146	Industrial Series Classic Die Cast Alloy Boss Head
04999147	Industrial Series Mini All Metal Boss Head Clamp
04999148	Industrial Series Heavy Duty Powder Coated Boss head
04999149	Industrial Series Universal Swivel Boss Head
04999150	Industrial Series Square Boss Had Clamp



© 2016 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at [www.fishersci.com/trademarks](http://www.fishersci.com/trademarks).

### In the United States:

For customer service, call 1-800-766-7000

To fax an order, use 1-800-926-1166

To order online: [www.fishersci.com](http://www.fishersci.com)

BN0116177

 **Fisher  
Scientific**

A Thermo Fisher Scientific Brand