

Blue

the choice is yours

White


By tradition, biopsy wrapping paper is white. When biopsy specimens are light in color this can result in an inability to clearly locate the specimen. If there are multiple fragments of tissue, there is a risk of losing or mislaying some of the pieces.

There needs to be a good contrast between the paper and the specimen, therefore we are pleased to add a blue TissueWrap paper to our range of Biopsy handling products.

CellPath TissueWrap Biopsy Paper is specifically designed for processing biopsy specimens. The tissue paper is cut to a convenient ready-to-use size.

The non-absorbent paper, available in white or blue, allows fluids to flow through freely while retaining small tissue fragments. TissueWrap is supplied in a dispenser pack for increased convenience.

- High chemical resistance
- Disposable – time saving
- Cost effective



For more information on any of the products in CellPath's range go to

www.cellpath.us.com

Call Toll Free on 1-855-CELPATH
Email us at sales@cellpath.us.com
Click Here to buy online

CellPath 
INNOVATION IN CELLULAR PATHOLOGY