



QUIDEL

Virena™



Global Wireless Surveillance and Remote Instrument Data Management



An innovative surveillance solution for the daily capture and management of patient de-identified results

An economical instrument management system that allows you the ability to collate, report and analyze daily test results from all your Sofia(s)

With Virena you can

- Generate dynamic reports, charts and trend assay usage
- Support operator training initiatives
- Monitor QC and calibration frequency
- Calculate percent positives
- Improve surveillance compliance year round and provide a community service by sharing de-identified surveillance results

*Sofia shown with Virena router.

How does Virena work?

Virena connects one or more Sofia(s) to a Virena router that transmits patient de-identified test results daily to a secure database. The Virena allows you the ability to compile each facility's results into a single database. In addition, Virena has an option for automatically sharing "near" real-time de-identified compliant results with public health. Participation in this program will be part of a growing national influenza surveillance program.



myvirena.com allows you to sort and group results

Trend data

Run charts and reports



QUIDEL

For more information or to order, call your Fisher HealthCare representative

© 2014 Quidel Corporation. All rights reserved. Quidel® and AmpliVue® trademarks used are owned by Quidel Corporation. Fisher HealthCare trademark used is owned as indicated at www.fishersci.com/trademarks.

BN1204141 | 1426M1014D (12/14)

In the United States:

For customer service, call 1-800-640-0640

To fax an order, use 1-800-290-0290

To order online: www.fisherhealthcare.com



Part of Thermo Fisher Scientific