

Supa Mega Adhesion Slides

NEW



HISTOBOND® + Supa Mega slides are designed for mounting large sections of neuro and breast tissue. Manufactured from the highest quality white glass, with an innovative positive charge coating, **HISTOBOND® + Supa Mega slides** offer exceptional adhesive properties and provide a superior coverage of reagents across the slide.

The HISTOBOND® + Supa Mega slides meet the needs of histopathology laboratories wanting to address problematic tissue specimens such as neuro, breast and fatty tissues, as well as orthopaedic specimens such as femoral heads. The excellent adhesive properties of the slide ensure that specimens remain securely attached during demanding antigen retrieval processes and special stain techniques.

The HISTOBOND® + Supa Mega slides have a 'Right Side' indicator (RS) to help ensure users correctly orientate the slide as well as a 7mm marking area across the width of the slide in order to maximise the working area.

- Approximate slide dimensions: 75.2 x 50.4 x 1mm
- Pack Size; 50
- Colours: White

Product Code	Unit Size	Description
MAD-1402-09A	50	HISTOBOND® + SUPA MEGA WHITE - MICROSCOPE SLIDE (CHARGED)



The brand Marienfeld Superior has become a synonym for quality and reliability, with experience of more than 89 years in developing and producing laboratory glassware.

For more information please contact your local Dealer



CellPath Ltd
80 Mochdre Enterprise Park
Newtown, Powys
SY16 4LE, Mid-Wales UK

Tel: +44 (0) 1686 611 333
Fax: +44 (0) 1686 622 946
Email: sales@cellpath.co.uk
Web: www.cellpath.co.uk