

# ProMarc

## Supa Mega Mother Ship

Accurate and reliable tracking of specimens is paramount to the histology laboratory. The ability to identify Supa Mega cassettes has, until now, been a process prone to human error.

CellPath's Supa Mega Mother Ship cassettes use an innovative permanent docking bay for a standard\* tissue processing cassette printed with a 2D bar code and/or test information.

- Significantly reduces turn-around-times.
- Improved specimen quality, reducing distortion with no curling or "pringle" effects.
- Improved processing and reduced carry-over due to a large open area of Hex pore design to maximize fluid exchange.

*\*Validated using a CellPath System II+ Processing Cassette*

00.849.930/1

**For more information please call 1-855-CELPATH**

**CellPath**   
INNOVATION IN CELLULAR PATHOLOGY

**CellPath USA**  
C/O Bonded Logistics  
5709 North Graham Street  
P.O. Box 480203  
Charlotte, NC 28269


Call: Toll Free 1-855-CELPATH  
Fax: Toll Free 1-855-CELPATH  
Email: [sales@cellpath.us.com](mailto:sales@cellpath.us.com)  
Web: [www.cellpath.us.com](http://www.cellpath.us.com)

# ProMarc

## Supa Mega Mother Ship



The CellTec Block Trimmer Plus allows users to remove excess paraffin wax efficiently from the outside of tissue cassettes, with no risk of damage to printed text or barcodes and with improved safety for laboratory staff. Melted paraffin wax is collected in a disposable collection pot.



The Block Trimmer Plus features an ExacTrim Locator designed to remove the wax from the writing surface of a standard cassette\* permanently docked in the Supa Mega Mother Ship. The ExacTrim can also be used as a fine detail trimmer for conventional blocks.

*The Supa Mega Mother Ship (patent pending) is available in 6 colors*

Product Code	Description	Pack Size
EAO-0102-02A	SUPA MEGA MOTHER SHIP - WHITE	100
EAO-0104-02A	SUPA MEGA MOTHER SHIP - GREEN	100
EAO-0106-02A	SUPA MEGA MOTHER SHIP - BLUE	100
EAO-0108-02A	SUPA MEGA MOTHER SHIP - PINK	100
EAO-0109-02A	SUPA MEGA MOTHER SHIP - YELLOW	100
EAO-0110-02A	SUPA MEGA MOTHER SHIP - RED	100
GBC-3614-05A	STAINLESS STEEL SUPA MEGA BASE MOULD 36x36x10mm	5
GBC-6014-05A	STAINLESS STEEL SUPA MEGA BASE MOULD 60x45x15mm	5
JAZ-0200-00A	BLOCK TRIMMER PLUS, NO TILT, 110 - 230V.a.c. 50-60Hz	1

**For more information please call 1-855-CELPATH**

**CellPath**   
 INNOVATION IN CELLULAR PATHOLOGY

**CellPath USA**  
 C/O Bonded Logistics  
 5709 North Graham Street  
 P.O. Box 480203  
 Charlotte, NC 28269

Call: Toll Free 1-855-CELPATH  
 Fax: Toll Free 1-855-CELPATH  
 Email: [sales@cellpath.us.com](mailto:sales@cellpath.us.com)  
 Web: [www.cellpath.us.com](http://www.cellpath.us.com)