

Discover. Synthesize. Analyze. Customize.

# New! Fisher BioReagents Granulated Microbiological Media

Introducing the most commonly used bacteriological growth media, now available in granulated form, which minimizes dusting and inhalation exposure. These media are carefully formulated for the propagation and maintenance of *E. coli* strains and are used in various molecular biology applications.

- Easy to weigh, dust-free granules
- Tested for growth promotion with various *E. coli* strains
- Shelf life is 3 years minimum when stored at room temperature

These products are now offered in three convenient pack sizes and polypropylene jars to meet your laboratory needs



Catalog No.	Description	Pack Size	Packaging
BP9743-500	2XTY Broth Granulated	500 g	Polypropylene jars
BP9743-2	2XTY Broth Granulated	2 kg	
BP9743-5	2XTY Broth Granulated	5 kg	
BP9744-500	Agar, Granulated	500 g	Polypropylene jars
BP9744-2	Agar, Granulated	2 kg	
BP9744-5	Agar, Granulated	5 kg	
BP9745-500	LB Agar, (Lennox L Agar), Granulated	500 g	Polypropylene jars
BP9745-2	LB Agar, (Lennox L Agar), Granulated	2 kg	

## For General Laboratory Use

- LB Agar (Lennox) formulation is suitable for the preparation of bacteriological pour plates
- 2XYT Broth is an enriched medium used for faster bacterial growth and for growth of hosts replicating M13 vectors
- Agar may be used as solidifying agent in bacteriological culture media

## Related Product

BP2820 Water, Microbial Cell Culture Grade

- 0.03 micron filtered and autoclaved to ensure purity
- Low level of endotoxins < 0.25 EU/mL
- Total organic carbon (TOC) < 50 ppb
- Functionally tested for promoting microbial growth



Fisher BioReagents™

© 2015 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at [www.fishersci.com/trademarks](http://www.fishersci.com/trademarks).

### In the United States:

For customer service, call 1-800-766-7000

To fax an order, use 1-800-926-1166

To order online: [www.fishersci.com](http://www.fishersci.com)

### In Canada:

For customer service, call 1-800-234-7437

To fax an order, use 1-800-463-2996

To order online: [www.fishersci.ca](http://www.fishersci.ca)



A Thermo Fisher Scientific Brand