

KIMAX® STARTER PACKS

• SAVE TIME • SAVE SPACE • SAVE MONEY WITH KIMAX® LAB STARTER PACKS

- Purchase only what you need
- Less storage space required
- The most popular sizes in one case saves time searching for the correct glassware



LAB STARTER PACKS



Graduated Cylinder Starter Pack, To Deliver

Catalog No.	Cap. (mL)	Grad. Interval (mL)	Subdiv. (mL)	±Tol. (mL)	Approx. Ht. (mm)	Case Qty.
03-340-197	10	1 to 10	0.2	0.1	135	1
	25	3 to 25	0.5	0.3	140	1
	50	3 to 50	1	0.4	190	1
	100	5 to 100	1	0.6	255	1
	250	10 to 250	2	1.4	330	1



Low-Form Beaker Starter Pack

Catalog No.	Cap. (mL)	Grad. Interval (mL)	Subdiv. (mL)	Approx. OD x Ht. (mm)	Case Qty.
02-555-2	50	10 to 40	10	42 x 53	1
	100	20 to 80	10	50 x 70	1
	250	25 to 200	10	68 x 88	1
	600	50 to 500	25	88 x 122	1
	1000	50 to 950	50	107 x 145	1



Heavy-Duty Beaker Starter Pack

Catalog No.	Cap. (mL)	Grad. Interval (mL)	Subdiv. (mL)	Approx. OD x Ht. (mm)	Case Qty.
02-555-101	250	25 to 200	25	67 x 92	1
	400	25 to 325	50	77 x 114	1
	600	50 to 500	50	87 x 122	1
	1000	100 to 1000	50	107 x 152	1



Erlenmeyer Flask Starter Pack

Catalog No.	Cap. (mL)	Grad. Interval (mL)	Subdiv. (mL)	Widest Dia. (mm)	Approx. Ht. (mm)	Case Qty.
10-102	50	20 to 50	10	50	78	1
	125	50 to 125	25	66	108	1
	250	50 to 225	25	81	130	1
	500	100 to 500	50	102	174	1
	1000	250 to 1000	50	128	213	1



Media Bottles Starter Pack

Catalog No.	Cap. (mL)	Grad. Interval (mL)	Subdiv. (mL)	Widest Dia. (mm)	Approx. Ht. (mm)	Case Qty.
02-542-355	50	20 to 50	10	50	78	1
	125	50 to 125	25	66	108	1
	250	50 to 225	25	81	130	1
	500	100 to 500	50	102	174	1
	1000	250 to 1000	50	128	213	1



© 2015 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at www.fishersci.com/trademarks.

In the United States:

For customer service, call 1-800-766-7000
 To fax an order, use 1-800-926-1166
 To order online: www.fishersci.com
 BN12091513

In Canada:

For customer service, call 1-800-234-7437
 To fax an order, use 1-800-463-2996
 To order online: www.fishersci.ca



A Thermo Fisher Scientific Brand