

The Altus™ table system provides just the right combination of flexibility, durability and affordability to meet current needs and allow for future changes. Altus™ is a functional, reliable and attractive laboratory workstation at competitive price.

# ALTUS™

## TABLE SYSTEM

**Highly Adaptable** - Flexible for changing laboratories; conveniently integrate these tables into your existing lab or use them as a basis of a freestanding easily modified new lab.

**Sustainable Labs** - Altus™ tables extend the productive life cycle of your laboratory by maximizing your ability to change locations, layouts and configurations.

**Improved Productivity** - The table frame can be quickly adjusted from 29" to 36" on 1" increments to easily adapt to new tasks, equipment and users.

**Durable** - Suitable for multiple applications and designed to withstand heavy loads.

**Adjustable Shelving** - Height adjustability in 1" increments delivers maximum adaptability to accommodate the requirements of individual users. Shelving options include powder coated steel, wood or glass.

**Materials, Coating Finish and Color** - Altus™ table systems are available in powder coated steel in a wide variety of standard or custom colors; stainless steel is an available option.

**Accessories** - Altus™ tables can be configured with many accessories such as keyboard trays, monitor arms, bulletin boards and shelf task lights to suit your needs.

**Storage Solutions** - Designed to accommodate suspended, mobile and floor mounted casework; reconfigure the storage space as your tasks change. Cabinets are available in steel, stainless steel, wood, and steel cabinets with wood fronts.



**Pre-Plumbing and Pre-Wiring** - Integrated plumbing, electrical and data allows for quick and easy access to services. Altus Table Systems are available UL Listed to UL 962 Standards.



**Overhead Service Solutions** - Altus™ tables are equipped to permit quick, easy, and affordable connection of electrical, mechanical, and data services for optimized layout flexibility. Mott offers many overhead service solutions to meet your project specific requirements.



**Cable Management Solutions** - Integral cable trays are a convenient and clutter reducing solution for power and data cable management.



**mott**  
MANUFACTURING

ISO 9001:2008 REGISTERED

QUALITY BY DESIGN

Corporate Headquarters - 452 Hardy Rd, Brantford ON, Canada N3T 5L8 • T. 519.752.7825  
Wood Casework Division - 562 Industrial Park Rd, P.O. Box 156, Maxwelton, WV, USA 24957 • T. 304.497.2115

[www.mott.ca](http://www.mott.ca)

January 2015



Altus shown with sliding glass wall cabinets and shelves  
Table frames positioned in sitting and standing height



Altus shown with bulletin boards



Mobile and suspended cabinets used for storage



**Mobile Solutions** - Heavy-duty swivel casters have an adjustment knob for times when you need both mobility and absolute stability.

**ALTUS™**  
TABLE SYSTEM