



Thermo Scientific Nalgene
Sterile Single-Use PETG Flasks



Culture with confidence



Thermo Scientific™ Nalgene™ Sterile Single-Use PETG Flasks are ideal for suspension cell culture, media preparation, mixing and storage.




These flasks are molded of BPA-free, crystal-clear, light-weight, break resistant PETG. Available with plain or baffled bottoms, vented or non-vented closures.

- Sterile and single use: Reduce contamination risk
- Robust and break resistant: Increase personal safety
- Collapses when autoclaved: Reduce biohazardous waste
- Graduated for quick volume assessment
- Excellent validation support: Validation binder available upon request*

* Validation binders are available to customers as support documentation for your product validation process. This confidential information can be requested at: rocrgsupport@thermofisher.com. Following your request submission, a Confidentiality Agreement will be forwarded to you for completion.

Ideal solution for scaling up your suspension cultures.

- Enhance your process with plug-and-play solutions
- Improve aseptic technique for your seed train
- Compatible with existing single-use process consumables



Nalgene Sterile PETG Single Use Erlenmeyer Flasks

Cat. No.	Working Capacity, mL	Bottom Style	Neck Finish, mm	No. Per Case
Solid closures				
10-041-20	125	Plain	38-430	24
10-346-7A	250	Plain	38-430	12
10-346-7B	500	Plain	45-430	12
10-041-21	1000	Plain	45-430	6
10-346-7C	2000	Plain	45-430	4
03-600-150	2800	Plain	70	4
10-041-23	125	Baffled	38-430	24
10-346-8A	250	Baffled	38-430	12
10-346-8B	500	Baffled	45-430	12
10-041-24	1000	Baffled	45-430	6
10-346-8C	2000	Baffled	45-430	4
03-600-151	2800	Baffled	70	4
Vented closures				
10-524	125	Plain	38-430	24
10-525	250	Plain	38-430	12
10-526	500	Plain	45-430	12
10-527	1000	Plain	45-430	6
10-528	2000	Plain	45-430	4
03-600-152	2800	Plain	70	4
10-529	125	Baffled	38-430	24
10-530	250	Baffled	38-430	12
10-531	500	Baffled	45-430	12
10-532	1000	Baffled	45-430	6
10-533	2000	Baffled	45-430	4
10-534	2800	Baffled	70	4
Replacement Vented Closures				
02-923-150	Vent Closures for 125 mL and 250 mL Erlenmeyer Flasks			12
02-923-151	Vent Closures for 500 mL, 1 L and 2 L Erlenmeyer Flasks			12

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.



In the United States:

For customer service, call 1-800-766-7000
To fax an order, use 1-800-926-1166
To order online: www.fishersci.com

BRRSDNALGSHKRFL 0516

In Canada:

For customer service, call 1-800-234-7437
To fax an order, use 1-800-463-2996
To order online: www.fishersci.com

BN0525164

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand