



Your One-Stop Shop

Food and beverage testing is essential. We understand the demands of this industry and want to provide you with the technologies and expertise needed to have confidence in your results.

Dispersive SPE: Efficient Sample Preparation and Cleanup

Thermo Scientific™ HyperSep™ Dispersive SPE Extraction Products are available in a range of pre-prepared sorbent combinations to help you save valuable preparation time, achieve consistent results, and more.



Extraction Tubes for Pesticide Residue Analysis Using QuEChers

Contents	Size	Quantity	Cat. No.
Magnesium Sulfate: 6g Sodium Acetate: 1.5g	50mL	250/Pack	03-150-631

Visit fishersci.com/hypersepextraction or fishersci.ca/hypersepextraction for more information.

Cleanup for the Schenck Variation of QuEChERS

Contents	Size	Quantity	Cat. No.
Carbon (200mg); PSA (400mg); Teflon Frit	6mL	30/Pack	60-105-207
Carbon (250mg); PSA (500mg); Teflon Frit	6mL	30/Pack	60-105-208
Carbon (500mg); PSA (500mg); Teflon Frit	6mL	30/Pack	60-105-209

Syringe Filters

Thermo Scientific™ Titan3™ and Target2™ syringe filters reliably eliminate both particles and microorganisms, protecting your columns and preventing premature blockages for more consistent experimental results.



Description	Diameter	Pore Size	Quantity	Cat. No.
Titan3 PTFE (Hydrophobic) Syringe Filters	30mm	1.0µm	100/Pack	14-823-282
Titan3 Nylon Syringe Filters	30mm	1.5µm	100/Pack	14-823-283
Target2 PVDF Syringe Filters	30mm	0.45µm	100/Pack	03-391-1F
Target2 Cellulose Acetate Syringe Filters	30mm	0.45µm	100/Pack	03-376-221

Note: Additional sizes and membrane types are available.

SPE Cartridges

Get rapid, effective, and targeted results from this broad range of versatile, high-quality Thermo Scientific™ HyperSep™ and HyperSep™ Hypercarb™ SPE cartridges and other popular products.



SPE Cartridges	Bed Weight	Column Capacity	Quantity	Cat. No.
HyperSep C18	200mg	3mL	50/Pack	03-251-259
HyperSep Hypercarb	50mg	1mL	50/Pack	03-250-400
HyperSep SCX	500mg	3mL	50/Pack	03-25-1275
HyperSep Silica	1g	6mL	50/Pack	03-25-1278

Note: Additional cartridge options are available.

Visit fishersci.com/thermochrom or fishersci.ca/thermochrom for more information.

Vials and Closures

Thermo Scientific™ SureStop™ Vials with Advanced Vial Closure System (AVCS) technology significantly reduce errors that result from evaporation of solvents in vials that are not sealed correctly.



Description	Closure Size	Volume	Quantity	Cat. No.
Autosampler Vial Kits, Unassembled	9mm	1.5mL	100/Pack	03-246-324
Certified Vial Kits	9mm	2.0mL	100/Pack	03-151-692

Note: Additional sizes available.

LC Columns

As a leader in HPLC column technology, Thermo Scientific products offer the highest resolution chemistries for your food testing analyses. Use the recommendation guide below to get started.



Test	Phase/ Chemistry	LC-UV/FLD/CAD/RI		LC-MS/MS	
		Size: 3µm, 3 x 150mm	Cat. No.	Size: 2.6µm, 2.1 x 100mm	Cat. No.
Mycotoxins	C18	Acclaim™ 120 C18 Columns	DX063691	Accucore™ aQ C18 Polar Endcapped	17-326-102130
Pesticides	Polar Embedded	Acclaim Polar Advantage II	DX063705		
Water Soluble	aQ	Hypersil™ GOLD aQ C18 Polar Endcapped	25-303-153030		
Fat Soluble Vitamins	C18	Hypersil GOLD C18 Selectivity	25-003-153030	Accucore RP-MS	17-626-102130
Lipids	C30	Acclaim C30	DX075724	Accucore C30	03-452-136
Sugars	HILIC	Hypersil GOLD Amino	25-703-103030	Accucore 150 Amide HILIC	03-452-147
Organic Acids	Mixed Mode	Acclaim LC Column	DX070086	Accucore RP-MS	17-626-102130

Visit fishersci.com/foodtesting or fishersci.ca/foodtesting to browse all food and beverage testing products.

In the United States

Order online: fishersci.com
 Fax an order: 1-800-926-1166
 Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca
 Fax an order: 1-800-463-2996
 Call customer service: 1-800-234-7437

