

# Honeywell Research Chemicals

## U.S. Portfolio Overview

Portfolio Brands

**Honeywell**  
Burdick & Jackson™

**Honeywell**  
Riedel-de Haën™

**Honeywell**  
**Fluka™**

**Honeywell**

Portfolio Overview

### SOLVENTS

High-purity solvents for multipurpose and dedicated applications

### SOLVENTS

High-purity solvents for dedicated applications

### ANALYTICAL REAGENTS

Karl Fischer titration, standards, acids, bases, salts, and pH buffers

### PERFORMANCE GRADE

Solvents and inorganics for general laboratory use

Product Lines

- **B&J Brand™** – multipurpose solvents for demanding analytical applications
- **B&J GC2™** – Capillary gas chromatography applications requiring trace analysis at or below the part-per-billion level
- **BioSyn™** – DNA, RNA and peptide synthesis
- **Anhydrous** – Moisture-sensitive organic synthesis, organometallic, combinatorial chemistry, and related applications
- **B&J Purified Plus™** – Specialty organic synthesis, prep-LC and other industrial and pharmaceutical applications requiring higher levels of lot-to-lot consistency

- **Chromasolv™** – Chromatography techniques including LC-MS, UHPLC, Headspace, and GC
- **TraceSELECT™ Solvents** – Trace and low-metal content analysis

- **Analytical Inorganics**
- **Standard Solutions** – Buffers and concentrated solutions
- **Hydranal™** – Karl Fischer titration for measuring water content
- **TraceSELECT™ Inorganics** – Trace and low-metal content analysis

- **Solvents** – ACS specific use (HPLC), ACS general use, solvents for histology, and reagent grade solvents for chemical synthesis and other industrial applications
- **Inorganics** – Chemical synthesis and inorganic chemistry including essential acids and bases, salts, metals and elements, and reagents for chemical reactions

Visit [fishersci.com/HoneywellChemicals](http://fishersci.com/HoneywellChemicals) to learn more, or visit [fishersci.com/ChemicalSpecialist](http://fishersci.com/ChemicalSpecialist) to connect with a Fisher Scientific Chemical Specialist who can answer your questions or help you place an order.

# Honeywell

© 2016 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at [www.fishersci.com/trademarks](http://www.fishersci.com/trademarks).

### In the United States:

For customer service, call 1-800-766-7000

To fax an order, use 1-800-926-1166

To order online: [fishersci.com](http://fishersci.com)

Burdick & Jackson, B&J Brand, B&J GC2, BioSyn, B&J Purified Plus, and Honeywell LabReady are trademarks of Honeywell International Inc.  
Riedel-de Haen, Fluka, Hydranal, Chromasolv, and TraceSELECT are trademarks of Honeywell Specialty Chemicals Seelze GmbH.

16-58482 AL 12/16 BN1209161



A Thermo Fisher Scientific Brand