



WHEATON® Glass Bottles

WHEATON offers a comprehensive line of glass laboratory bottles and jars. Glass bottles offer sample reliability and integrity, and is a mainstay of labs worldwide. WHEATON bottles are manufactured from USP Type III soda-lime glass with moderate chemical resistance. Popular bottle styles available include Media Bottles, Boston Rounds, Wide Mouth Packers, Straight Sided Jars and Safety Coated Bottles for enhanced laboratory protection. WHEATON completes the package with the right closure. We provide a variety of caps and seals to ensure a perfect fit for each container.

WHEATON offers convenience bulk packs of containers with or without caps attached for high use items or facilities with centralized stockrooms. With safety in mind, Tamper Evident Seal / HAZCOM Labels are provided in each case of bottles with the exception of bulk packs.

Product Highlights	
<ul style="list-style-type: none"> Bottles manufactured from clear and amber USP Type III soda-lime glass 	<ul style="list-style-type: none"> Wide mouth bottles are ideal for dry and viscous samples
<ul style="list-style-type: none"> Amber colored glass is ideal for light sensitive products 	<ul style="list-style-type: none"> White polypropylene or black phenolic screw closures can be purchased separately
<ul style="list-style-type: none"> Narrow mouth bottles are ideal for liquids 	<ul style="list-style-type: none"> Methods of sterilization: dry heat or EtO (not suitable for autoclave)
<ul style="list-style-type: none"> Safety coated glass prevents dangerous spills and control sharp fragments 	

French Square

- Clear or Amber, USP Type III soda-lime glass
- Ideal for solvent, chemical or sample storage
- Square shape maximizes storage space
- Available with caps attached or bulk packs without caps

Valumetric™ Graduated Bottle

- Clear, USP Type III soda-lime glass
- Ideal for general storage
- Easy removal of contents
- Graduated in mL and ounces
- Available with caps attached or bulk packs without caps

Media Bottle, Lab 45

- Manufactured from USP Type I borosilicate glass
- Available with or without caps
- 45mm screw thread, wide mouth finish
- Cap has a built-in wedge shaped sealing, eliminating the need for a liner



Boston Round Bottle

- Clear or Amber, Type III soda-lime glass
- Ideal for solvent, chemical or sample storage
- Available with caps attached or bulk packs without caps

Media Bottles

- Clear borosilicate glass bottle with graduations and marking patch
- Amber borosilicate glass bottle in 125mL size only
- Pre-attached black phenolic screw caps with choice of three liners
- Hand grips on 1 Liter bottle provide easy handling

Wide Mouth Packer

- Amber, USP Type III soda-lime glass
- Ideal for storage of light sensitive samples or dry powders
- Easy removal of contents

© 2014 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at www.fishersci.com/trademarks.

In the United States:

For customer service, call 1-800-766-7000
To fax an order, use 1-800-926-1166
To order online: www.fishersci.com

In Canada:

For customer service, call 1-800-234-7437
To fax an order, use 1-800-463-2996
To order online: www.fishersci.ca



A Thermo Fisher Scientific Brand