

Your guide to chromatography

Pure proteins today. Powerful results tomorrow

Click on the chromatography method you are interested in to jump to the corresponding selection guide.







AC
Affinity chromatography

SEC
Size exclusion chromatography

IEX
Ion exchange chromatography




GE format guide

Type of purification	Manual purification			Manual or system purification		System purification
Format	 Spin columns	 96-well plates	 Gravity flow columns	 Bottles of chromatography media (resins)	 Small-column cartridges	 Other columns
Format name	SpinTrap™	MultiTrap™	GraviTrap™, MiniTrap™, MidiTrap™, PD10	Lab pack	HiTrap™	HiScreen™, HiPrep™, HiLoad™, RESOURCE™, Tricorn™, and Precision
Use	Screening and quick desalting of small sample quantities using a benchtop centrifuge	High-throughput screening and small scale purification using centrifuge or vacuum equipment	Simple one-step purification of proteins or sample desalting without the need for equipment	Batch purification and self-packing	Easy to use with a syringe, peristaltic pump, or a chromatography system	Larger scale or high performance applications

Affinity chromatography

Guide to chromatography media


Click on the icons to be directed to the corresponding web page.



Size exclusion chromatography

Guide to chromatography media

Click on the icons to be directed to the corresponding web page.




Ion exchange chromatography

Home

Guide to chromatography media

Click on the icons to be directed to the corresponding web page.



Need further guidance?

Download our free Purify App

<https://www.fishersci.com/us/en/brands/I9C8LWZ9/resources.html>



For ordering information, visit www.fishersci.com/gehc

GE, GE monogram, Cpto, GraviTrap, HiLoad, HiPrep, HiScreen, HiTrap, MabSelect SuRe, MidiTrap, MiniBeads, MiniTrap, MonoBeads, MultiTrap, RESOURCE, Sephacryl, Sephadex, Sepharose, SOURCE, SpinTrap, Superdex, Superose, and Tricorn are trademarks of General Electric Company. Strep-tag and Strep-tactin are trademarks of IBA GmbH. TALON is a trademark of Clontech Laboratories, Inc. Superflow is a trademark of Sterogene Bioseparations Inc. All other trademarks are the property of their respective owners.
© 2015 General Electric Company. First published Aug. 2015.
GE Healthcare UK Limited, Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK.



Distributor
GE Healthcare



A Thermo Fisher Scientific Brand

© 2015 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at www.fishersci.com/trademarks.

In the United States:

For customer service, call 1-800-766-7000
To fax an order, use 1-800-926-1166
To order online: www.fishersci.com

In Canada:

For customer service, call 1-800-234-7437
To fax an order, use 1-800-463-2996
To order online: www.fishersci.ca