

Acrylic Lab Workspace Partitions

- Help minimize workspace proximity issues in your laboratory
- These clear, acrylic partitions allow the passage of light to the work area, providing for a safe, bright work environment
- Partitions are easy to clean and disinfect
- Our workspace partitions are available in 2 standard designs:
 - **Perimeter:** to laterally separate a workspace
 - **Island:** to provide a 3-sided separation barrier

Perimeter Partition

Cat. No.: 06000908



- All partitions are constructed from 1/4" acrylic
- All partitions include lower cord/hose passage points
- Perimeter style partitions are 24" tall x 24" deep and have triangular stabilizing feet

Island Partition

Cat. No.: 06000909



- Island style partitions are 24" tall x 24" deep x 36" wide and include stabilization feet and rear acrylic cleats for a secure assembly
- All units require some assembly

Custom sizes and shapes can be quickly designed & built to suit your needs

Contact Your Fisher Sales Representative for Support

Find out more at fishersci.com/fisherbrandlabfurnishings

In the United States

Order online: fishersci.com
Fax an order: 1-800-926-1166
Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca
Fax an order: 1-800-463-2996
Call customer service: 1-800-234-7437