

# Fisherbrand Locking Test Tube Rack

## Empty Tube Contents Quickly and Easily

The Fisherbrand™ Locking Test Tube Rack features a silicon grip that securely holds tubes in place while pouring.

- Holds 13mm and 16mm test tubes
- Securely grips up to 50 tubes
- Made from blue acrylonitrile butadiene styrene



### Ordering Information

Product	Cat. No.
Locking Test Tube Rack 50 Positions (10 x 5)	14-955-035



© 2017 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at [fishersci.com/trademarks](http://fishersci.com/trademarks).

#### In the United States:

For customer service, call 1-800-766-7000  
To fax an order, use 1-800-926-1166  
To order online: [fishersci.com](http://fishersci.com)

#### In Canada:

For customer service, call 1-800-234-7437  
To fax an order, use 1-800-463-2996  
To order online: [fishersci.ca](http://fishersci.ca)