

Standard Methods for the Examination of Water and Wastewater
22nd Edition

TABLE OF CONTENTS

Part 1000 INTRODUCTION

- 1010 INTRODUCTION
 - A. Scope and Application of Methods
 - B. Statistics
 - C. Glossary
 - D. Dilution/Concentration Options
- 1020 QUALITY ASSURANCE
 - A. Introduction
 - B. Quality Control
 - C. Quality Assessment
- 1030 DATA QUALITY
 - A. Introduction
 - B. Measurement Uncertainty
 - C. Method Detection Level
 - D. Data Quality Objectives
 - E. Checking Analyses' Correctness
- 1040 METHOD DEVELOPMENT AND EVALUATION
 - A. Introduction
 - B. Method Validation
 - C. Collaborative Testing
- 1050 EXPRESSION OF RESULTS
 - A. Units
 - B. Significant Figures
- 1060 COLLECTION AND PRESERVATION OF SAMPLES
 - A. Introduction
 - B. Collection of Samples
 - C. Sample Storage and Preservation

- 1080 REAGENT WATER
 - A. Introduction
 - B. Methods for Preparing Reagent-Grade Water
 - C. Reagent Water Quality
- 1090 LABORATORY OCCUPATIONAL HEALTH AND SAFETY
 - A. Introduction
 - B. Safe Laboratory Practices
 - C. Laboratory Facility/Fixed Equipment
 - D. Hazard Evaluation
 - E. Personal Protective Equipment
 - F. Worker Protection Medical Program
 - G. Provisions for Work with Particularly Hazardous Substances
 - H. Biological Safety
 - I. Radiological Safety
 - J. Chemical Hygiene Plan
- 1100 WASTE MINIMIZATION AND DISPOSAL
 - A. Introduction
 - B. Waste Minimization
 - C. Waste Treatment and Disposal

Part 2000 PHYSICAL AND AGGREGATE PROPERTIES

- 2010 INTRODUCTION
- 2020 QUALITY ASSURANCE/QUALITY CONTROL
 - A. Introduction
 - B. Quality Control Practices
- 2110 APPEARANCE
- 2120 COLOR
 - A. Introduction
 - B. Visual Comparison Method
 - C. Spectrophotometric—Single-Wavelength Method (Proposed)
 - D. Spectrophotometric—Multi-Wavelength Method
 - E. Tristimulus Spectrophotometric Method
 - F. ADMI Weighted-Ordinate Spectrophotometric Method
- 2130 TURBIDITY
 - A. Introduction
 - B. Nephelometric Method
- 2150 ODOR
 - A. Introduction
 - B. Threshold Odor Test
- 2160 TASTE
 - A. Introduction
 - B. Flavor Threshold Test (FTT)
 - C. Flavor Rating Assessment (FRA)
- 2170 FLAVOR PROFILE ANALYSIS
 - A. Introduction
 - B. Flavor Profile Analysis

- 2310 ACIDITY
 - A. Introduction
 - B. Titration Method
- 2320 ALKALINITY
 - A. Introduction
 - B. Titration Method
- 2330 CALCIUM CARBONATE SATURATION
 - A. Introduction
 - B. Indices Indicating Tendency of a Water to Precipitate or Dissolve CaCO_3
 - C. Indices Predicting the Quantity of CaCO_3 That Can Be Precipitated or Dissolved
 - D. Graphical and Computer Methods for CaCO_3 Indices
- 2340 HARDNESS
 - A. Introduction
 - B. Hardness by Calculation
 - C. EDTA Titrimetric Method
- 2350 OXIDANT DEMAND/REQUIREMENT
 - A. Introduction
 - B. Chlorine Demand/Requirement
 - C. Chlorine Dioxide Demand/Requirement
 - D. Ozone Demand/Requirement—Batch Method
 - E. Ozone Demand/Requirement—Semi-Batch Method
- 2510 CONDUCTIVITY
 - A. Introduction
 - B. Laboratory Method
- 2520 SALINITY
 - A. Introduction
 - B. Electrical Conductivity Method
 - C. Density Method
 - D. Algorithm of Practical Salinity
- 2530 FLOATABLES
 - A. Introduction
 - B. Particulate Floatables
 - C. Trichlorotrifluoroethane-Soluble Floatable Oil and Grease
- 2540 SOLIDS
 - A. Introduction
 - B. Total Solids Dried at 103–105°C
 - C. Total Dissolved Solids Dried at 180°C
 - D. Total Suspended Solids Dried at 103–105°C
 - E. Fixed and Volatile Solids Ignited at 550°C
 - F. Settleable Solids
 - G. Total, Fixed, and Volatile Solids in Solid and Semisolid Samples
- 2550 TEMPERATURE
 - A. Introduction
 - B. Laboratory and Field Methods
- 2560 PARTICLE COUNTING AND SIZE DISTRIBUTION
 - A. Introduction

- B. Electrical Sensing Zone Method
 - C. Light-Blockage Methods
 - D. Light-Scattering Method
- 2570 ASBESTOS
- A. Introduction
 - B. Transmission Electron Microscopy Method
- 2580 OXIDATION–REDUCTION POTENTIAL (ORP)
- A. Introduction
 - B. Oxidation–Reduction Potential Measurement in Clean Water
- 2710 TESTS ON SLUDGES
- A. Introduction
 - B. Oxygen-Consumption Rate
 - C. Settled Sludge Volume
 - D. Sludge Volume Index
 - E. Zone Settling Rate
 - F. Specific Gravity
 - G. Capillary Suction Time
 - H. Time-to-Filter
 - I. Modified Settled Sludge Volume
- 2720 ANAEROBIC SLUDGE DIGESTER GAS ANALYSIS
- A. Introduction
 - B. Volumetric Method
 - C. Gas Chromatographic Method
- 2810 DISSOLVED GAS SUPERSATURATION
- A. Introduction
 - B. Direct-Sensing Membrane-Diffusion Method

Part 3000 METALS

- 3010 INTRODUCTION
- A. General Discussion
 - B. Sampling and Sample Preservation
 - C. General Precautions
- 3020 QUALITY ASSURANCE/QUALITY CONTROL
- A. Introduction
 - B. Quality Control Practices
- 3030 PRELIMINARY TREATMENT OF SAMPLES
- A. Introduction
 - B. Filtration for Dissolved and Suspended Metals
 - C. Treatment for Acid-Extractable Metals
 - D. Digestion for Metals
 - E. Nitric Acid Digestion
 - F. Nitric Acid-Hydrochloric Acid Digestion
 - G. Nitric Acid-Sulfuric Acid Digestion
 - H. Nitric Acid-Perchloric Acid Digestion
 - I. Nitric Acid-Perchloric Acid-Hydrofluoric Acid Digestion
 - J. Dry Ashing

- K. Microwave-Assisted Digestion
- 3110 METALS BY ATOMIC ABSORPTION SPECTROMETRY
- 3111 METALS BY FLAME ATOMIC ABSORPTION SPECTROMETRY
 - A. Introduction
 - B. Direct Air-Acetylene Flame Method
 - C. Extraction/Air-Acetylene Flame Method
 - D. Direct Nitrous Oxide-Acetylene Flame Method
 - E. Extraction/Nitrous Oxide-Acetylene Flame Method
- 3112 METALS BY COLD-VAPOR ATOMIC ABSORPTION SPECTROMETRY
 - A. Introduction
 - B. Cold-Vapor Atomic Absorption Spectrometric Method
- 3113 METALS BY ELECTROTHERMAL ATOMIC ABSORPTION SPECTROMETRY
 - A. Introduction
 - B. Electrothermal Atomic Absorption Spectrometric Method
- 3114 ARSENIC AND SELENIUM BY HYDRIDE GENERATION/ATOMIC ABSORPTION SPECTROMETRY
 - A. Introduction
 - B. Manual Hydride Generation/Atomic Absorption Spectrometric Method
 - C. Continuous Hydride Generation/Atomic Absorption Spectrometric Method
- 3120 METALS BY PLASMA EMISSION SPECTROSCOPY
 - A. Introduction
 - B. Inductively Coupled Plasma (ICP) Method
- 3125 METALS BY INDUCTIVELY COUPLED PLASMA/MASS SPECTROMETRY
 - A. Introduction
 - B. Inductively-Coupled Plasma/Mass Spectrometry (ICP/MS) Method
- 3130 METALS BY ANODIC STRIPPING VOLTAMMETRY
 - A. Introduction
 - B. Determination of Lead, Cadmium, and Zinc
- 3500-Al ALUMINUM
 - A. Introduction
 - B. Eriochrome Cyanine R Method
- 3500-As ARSENIC
 - A. Introduction
 - B. Silver Diethyldithiocarbamate Method
- 3500-Ca CALCIUM
 - A. Introduction
 - B. EDTA Titrimetric Method
- 3500-Cr CHROMIUM
 - A. Introduction
 - B. Colorimetric Method
 - C. Ion Chromatographic Method
- 3500-Cu COPPER
 - A. Introduction
 - B. Neocuproine Method
 - C. Bathocuproine Method
- 3500-Fe IRON
 - A. Introduction

- B. Phenanthroline Method
- 3500-Pb LEAD
 - A. Introduction
 - B. Dithizone Method
- 3500-Li LITHIUM
 - A. Introduction
 - B. Flame Emission Photometric Method
- 3500-Mg MAGNESIUM
 - A. Introduction
 - B. Calculation Method
- 3500-Mn MANGANESE
 - A. Introduction
 - B. Persulfate Method
- 3500-K POTASSIUM
 - A. Introduction
 - B. Flame Photometric Method
 - C. Potassium-Selective Electrode Method
- 3500-Se SELENIUM
 - A. Introduction
 - B. Sample Preparation
 - C. Colorimetric Method
 - D. Determination of Volatile Selenium
 - E. Determination of Nonvolatile Organic Selenium Compounds
- 3500-Na SODIUM
 - A. Introduction
 - B. Flame Emission Photometric Method
- 3500-Sr STRONTIUM
 - A. Introduction
 - B. Flame Emission Photometric Method
- 3500-V VANADIUM
 - A. Introduction
 - B. Gallic Acid Method
- 3500-Zn ZINC
 - A. Introduction
 - B. Zincon Method
- 3500 OTHER METALS
 - 3500-Sb Antimony
 - 3500-Ba Barium
 - 3500-Be Beryllium
 - 3500-Bi Bismuth
 - 3500-B Boron
 - 3500-Cd Cadmium
 - 3500-Cs Cesium
 - 3500-Co Cobalt
 - 3500-Ga Gallium
 - 3500-Ge Germanium

3500-Au	Gold
3500-In	Indium
3500-Ir	Iridium
3500-Hg	Mercury
3500-Mo	Molybdenum
3500-Ni	Nickel
3500-Os	Osmium
3500-Pd	Palladium
3500-Pt	Platinum
3500-Re	Rhenium
3500-Rh	Rhodium
3500-Ru	Ruthenium
3500-Ag	Silver
3500-Te	Tellurium
3500-Tl	Thallium
3500-Th	Thorium
3500-Sn	Tin
3500-Ti	Titanium
3500-U	Uranium

Part 4000 INORGANIC NONMETALLIC CONSTITUENTS

4010 INTRODUCTION

4020 QUALITY ASSURANCE/QUALITY CONTROL

- A. Introduction
- B. Quality Control Practices

4110 DETERMINATION OF ANIONS BY ION CHROMATOGRAPHY

- A. Introduction
- B. Ion Chromatography with Chemical Suppression of Eluent Conductivity
- C. Single-Column Ion Chromatography with Direct Conductivity Detection
- D. Ion Chromatographic Determination of Oxyhalides and Bromide

4120 SEGMENTED CONTINUOUS FLOW ANALYSIS

- A. Introduction
- B. Segmented Flow Analysis Method

4130 INORGANIC NONMETALS BY FLOW INJECTION ANALYSIS

- A. Introduction
- B. Quality Control

4140 INORGANIC ANIONS BY CAPILLARY ION ELECTROPHORESIS

- A. Introduction
- B. Capillary Ion Electrophoresis with Indirect UV Detection

4500-B BORON

- A. Introduction
- B. Curcumin Method
- C. Carmine Method

4500-Br⁻ BROMIDE

- A. Introduction
- B. Phenol Red Colorimetric Method

- C. (Reserved)
- D. Flow Injection Analysis

4500-CO₂

CARBON DIOXIDE

- A. Introduction
- B. Nomographic Determination of Free Carbon Dioxide and the Three Forms of Alkalinity
- C. Titrimetric Method for Free Carbon Dioxide
- D. Carbon Dioxide and Forms of Alkalinity by Calculation

4500-CN⁻

CYANIDE

- A. Introduction
- B. Preliminary Treatment of Samples
- C. Total Cyanide after Distillation
- D. Titrimetric Method
- E. Colorimetric Method
- F. Cyanide-Selective Electrode Method
- G. Cyanides Amenable to Chlorination after Distillation
- H. Cyanides Amenable to Chlorination without Distillation (Short-Cut Method)
 - I. Weak Acid Dissociable Chloride
 - J. Cyanogen Chloride
- K. Spot Testing for Sample Screening
- L. Cyanates
- M. Thiocyanate
- N. Total Cyanide after Distillation, by Flow Injection Analysis
- O. Total Cyanide and Weak Acid Dissociable Cyanide by Flow Injection Analysis

4500-Cl CHLORINE (RESIDUAL)

- A. Introduction
- B. Iodometric Method I
- C. Iodometric Method II
- D. Amperometric Titration Method
- E. Low-Level Amperometric Titration Method
- F. DPD Ferrous Titrimetric Method
- G. DPD Colorimetric Method
- H. Syringaldazine (FACTS) Method
- I. Iodometric Electrode Technique

4500-Cl⁻

CHLORIDE

- A. Introduction
- B. Argentometric Method
- C. Mercuric Nitrate Method
- D. Potentiometric Method
- E. Automated Ferricyanide Method
- F. (Reserved)
- G. Mercuric Thiocyanate Flow Injection Analysis

4500-ClO₂

CHLORINE DIOXIDE

- A. Introduction
- B. Iodometric Method

- C. Amperometric Method I
 - D. (Reserved)
 - E. Amperometric Method II
- 4500-F⁻ FLUORIDE
- A. Introduction
 - B. Preliminary Distillation Step
 - C. Ion-Selective Electrode Method
 - D. SPADNS Method
 - E. Complexone Method
 - F. (Reserved)
 - G. Ion-Selective Electrode Flow Injection Analysis
- 4500-H⁺ PH VALUE
- A. Introduction
 - B. Electrometric Method
- 4500-I IODINE
- A. Introduction
 - B. Leuco Crystal Violet Method
 - C. Amperometric Titration Method
- 4500-I⁻ IODIDE
- A. Introduction
 - B. Leuco Crystal Violet Method
 - C. Catalytic Reduction Method
 - D. Voltammetric Method
- 4500-IO₃⁻ IODATE
- A. Introduction
 - B. Polarographic Method
- 4500-N NITROGEN
- A. Introduction
 - B. In-Line UV/Persulfate Digestion and Oxidation with Flow Injection Analysis
 - C. Persulfate Method
 - D. Conductimetric Determination of Inorganic Nitrogen
- 4500-NH₃ NITROGEN (AMMONIA)
- A. Introduction
 - B. Preliminary Distillation Step
 - C. Titrimetric Method
 - D. Ammonia-Selective Electrode Method
 - E. Ammonia-Selective Electrode Method Using Known Addition
 - F. Phenate Method
 - G. Automated Phenate Method
 - H. Flow Injection Analysis
- 4500-NO₂⁻ NITROGEN (NITRITE)
- A. Introduction
 - B. Colorimetric Method
- 4500-NO₃⁻ NITROGEN (NITRATE)
- A. Introduction
 - B. Ultraviolet Spectrophotometric Screening Method

- C. Second-Derivative Ultraviolet Spectrophotometric Method (Proposed)
 - D. Nitrate Electrode Method
 - E. Cadmium Reduction Method
 - F. Automated Cadmium Reduction Method
 - G. (Reserved)
 - H. Automated Hydrazine Reduction Method
 - I. Cadmium Reduction Flow Injection Analysis
- 4500-N_{org} NITROGEN (ORGANIC)
- A. Introduction
 - B. Macro-Kjeldahl Method
 - C. Semi-Micro Kjeldahl Method
 - D. Block Digestion and Flow Injection Analysis
- 4500-O OXYGEN (DISSOLVED)
- A. Introduction
 - B. Iodometric Methods
 - C. Azide Modification
 - D. Permanganate Modification
 - E. Alum Flocculation Modification
 - F. Copper Sulfate-Sulfamic Acid Flocculation Modification
 - G. Membrane Electrode Method
- 4500-O₃ OZONE (RESIDUAL)
- A. Introduction
 - B. Indigo Colorimetric Method
- 4500-P PHOSPHORUS
- A. Introduction
 - B. Sample Preparation
 - C. Vanadomolybdophosphoric Acid Colorimetric Method
 - D. Stannous Chloride Method
 - E. Ascorbic Acid Method
 - F. Automated Ascorbic Acid Reduction Method
 - G. Flow Injection Analysis for Orthophosphate
 - H. Manual Digestion and Flow Injection Analysis for Total Phosphorus
 - I. In-line UV/Persulfate Digestion and Flow Injection Analysis for Total Phosphorus
 - J. Persulfate Method for Simultaneous Determination of Total Nitrogen and Total Phosphorus
- 4500-KMnO₄ POTASSIUM PERMANGANATE
- A. Introduction
 - B. Spectrophotometric Method
- 4500-SiO₂ SILICA
- A. Introduction
 - B. (Reserved)
 - C. Molybdosilicate Method
 - D. Heteropoly Blue Method
 - E. Automated Method for Molybdate-Reactive Silica
 - F. Flow Injection Analysis for Molybdate-Reactive Silica
- 4500-S²⁻ SULFIDE
- A. Introduction

- B. Separation of Soluble and Insoluble Sulfides
- C. Sample Pretreatment to Remove Interfering Substances or to Concentrate the Sulfide
- D. Methylene Blue Method
- E. Gas Dialysis, Automated Methylene Blue Method
- F. Iodometric Method
- G. Ion-Selective Electrode Method
- H. Calculation of Un-ionized Hydrogen Sulfide
- I. Distillation, Methylene Blue Flow Injection Analysis
- J. Acid-Volatile Sulfide

4500-SO₃²⁻

SULFITE

- A. Introduction
- B. Iodometric Method
- C. Phenanthroline Method

4500-SO₄²⁻

SULFATE

- A. Introduction
- B. (Reserved)
- C. Gravimetric Method with Ignition of Residue
- D. Gravimetric Method with Drying of Residue
- E. Turbidimetric Method
- F. Automated Methylthymol Blue Method
- G. Methylthymol Blue Flow Injection Analysis

Part 5000 AGGREGATE ORGANIC CONSTITUENTS

5010 INTRODUCTION

- A. General Discussion
- B. Sample Collection and Preservation

5020 QUALITY ASSURANCE/QUALITY CONTROL

- A. Introduction
- B. Quality Control Practices

5210 BIOCHEMICAL OXYGEN DEMAND (BOD)

- A. Introduction
- B. 5-Day BOD Test
- C. Ultimate BOD Test
- D. Respirometric Method

5220 CHEMICAL OXYGEN DEMAND (COD)

- A. Introduction
- B. Open Reflux Method
- C. Closed Reflux, Titrimetric Method
- D. Closed Reflux, Colorimetric Method

5310 TOTAL ORGANIC CARBON (TOC)

- A. Introduction
- B. High-Temperature Combustion Method
- C. Persulfate-Ultraviolet or Heated-Persulfate Oxidation Method
- D. Wet Oxidation Method

5320 DISSOLVED ORGANIC HALOGEN

- A. Introduction

- B. Adsorption-Pyrolysis-Titrimetric Method
- 5510 AQUATIC HUMIC SUBSTANCES
 - A. Introduction
 - B. Diethylaminoethyl (DEAE) Method
 - C. XAD Method
- 5520 OIL AND GREASE
 - A. Introduction
 - B. Liquid-Liquid, Partition-Gravimetric Method
 - C. Partition-Infrared Method
 - D. Soxhlet Extraction Method
 - E. Extraction Method for Sludge Samples
 - F. Hydrocarbons
 - G. Solid-Phase, Partition-Gravimetric Method
- 5530 PHENOLS
 - A. Introduction
 - B. Cleanup Procedure
 - C. Chloroform Extraction Method
 - D. Direct Photometric Method
- 5540 SURFACTANTS
 - A. Introduction
 - B. Surfactant Separation by Sublation
 - C. Anionic Surfactants as MBAS
 - D. Nonionic Surfactants as CTAS
- 5550 TANNIN AND LIGNIN
 - A. Introduction
 - B. Colorimetric Method
- 5560 ORGANIC AND VOLATILE ACIDS
 - A. Introduction
 - B. Chromatographic Separation Method for Organic Acids
 - C. Distillation Method
 - D. Gas Chromatographic Method
- 5710 FORMATION OF TRIHALOMETHANES AND OTHER DISINFECTION BY-PRODUCTS
 - A. Introduction
 - B. Trihalomethane Formation Potential (THMFP)
 - C. Simulated Distribution System Trihalomethanes (SDS-THM)
 - D. Formation of Other Disinfection By-Products (DBPs)
- 5910 UV-ABSORBING ORGANIC CONSTITUENTS
 - A. Introduction
 - B. Ultraviolet Absorption Method

Part 6000 INDIVIDUAL ORGANIC COMPOUNDS

- 6010 INTRODUCTION
 - A. General Discussion
 - B. Sample Collection and Preservation
 - C. Analytical Methods
- 6020 QUALITY ASSURANCE/QUALITY CONTROL

- A. Introduction
- B. Quality Control Practices
- 6040 CONSTITUENT CONCENTRATION BY GAS EXTRACTION
 - A. Introduction
 - B. Closed-Loop Stripping, Gas Chromatographic/Mass Spectrometric Analysis
 - C. Purge and Trap Technique
 - D. Solid-Phase Microextraction (SPME)
 - E. Solid-Phase Microextraction (SPME) with CI GC/MS/MS
- 6200 VOLATILE ORGANIC COMPOUNDS
 - A. Introduction
 - B. Purge and Trap Capillary-Column Gas Chromatographic/Mass Spectrometric Method
 - C. Purge and Trap Capillary-Column Gas Chromatographic Method
- 6211 METHANE
 - A. Introduction
 - B. Combustible-Gas Indicator Method
 - C. Volumetric Method
- 6231 1,2-DIBROMOETHANE (EDB) AND 1,2-DIBROMO-3-CHLOROPROPANE (DBCP)
 - A. Introduction
 - B. Liquid-Liquid Extraction Gas Chromatographic Method
 - C. Purge and Trap Gas Chromatographic/Mass Spectrometric Method
 - D. Purge and Trap Gas Chromatographic Method
- 6232 TRIHALOMETHANES AND CHLORINATED ORGANIC SOLVENTS
 - A. Introduction
 - B. Liquid-Liquid Extraction Gas Chromatographic Method
 - C. Purge and Trap Gas Chromatographic/Mass Spectrometric Method
 - D. Purge and Trap Gas Chromatographic Method
- 6251 DISINFECTION BY-PRODUCTS: HALOACETIC ACIDS AND TRICHLOROPHENOL
 - A. Introduction
 - B. Micro Liquid-Liquid Extraction Gas Chromatographic Method
- 6252 DISINFECTION BY-PRODUCTS: ALDEHYDES
 - A. Introduction
 - B. PFBHA Liquid-Liquid Extraction Gas Chromatographic Method
- 6410 EXTRACTABLE BASE/NEUTRALS AND ACIDS
 - A. Introduction
 - B. Liquid-Liquid Extraction Gas Chromatographic/Mass Spectrometric Method
- 6420 PHENOLS
 - A. Introduction
 - B. Liquid-Liquid Extraction Gas Chromatographic Method
 - C. Liquid-Liquid Extraction Gas Chromatographic/Mass Spectrometric Method
- 6431 POLYCHLORINATED BIPHENYLS (PCBs)
 - A. Introduction
 - B. Liquid-Liquid Extraction Gas Chromatographic Method
 - C. Liquid-Liquid Extraction Gas Chromatographic/Mass Spectrometric Method
- 6440 POLYNUCLEAR AROMATIC HYDROCARBONS
 - A. Introduction
 - B. Liquid-Liquid Extraction Chromatographic Method

- C. Liquid-Liquid Extraction Gas Chromatographic/Mass Spectrometric Method
- 6450 NITROSAMINES
 - A. Introduction
 - B. Carbonaceous-Resin Solid-Phase Extraction GC/MS Method
 - C. Micro Liquid-Liquid Extraction GC/MS Method
- 6610 CARBAMATE PESTICIDES
 - A. Introduction
 - B. High-Performance Liquid Chromatographic Method
- 6630 ORGANOCHLORINE PESTICIDES
 - A. Introduction
 - B. Liquid-Liquid Extraction Gas Chromatographic Method I
Appendix—Standardization of Magnesia-Silica Gel Column by Weight Adjustment Based on Adsorption of Lauric Acid
 - C. Liquid-Liquid Extraction Gas Chromatographic Method II
 - D. Liquid-Liquid Extraction Gas Chromatographic/Mass Spectrometric Method
- 6640 ACIDIC HERBICIDE COMPOUNDS
 - A. Introduction
 - B. Micro Liquid-Liquid Extraction Gas Chromatographic Method
- 6651 GLYPHOSATE HERBICIDE
 - A. Introduction
 - B. Liquid Chromatographic Post-Column Fluorescence Method
- 6710 TRIBUTYL TIN
 - A. Introduction
 - B. Gas Chromatographic/Mass Spectrometric Method
 - C. Gas Chromatographic/Flame Photometric Detector Method

Part 7000 RADIOACTIVITY

- 7010 INTRODUCTION
 - A. General Discussion
 - B. Sample Collection and Preservation
- 7020 QUALITY SYSTEM
 - A. Quality Systems/Quality Assurance/Quality Control Program
 - B. Quality Control for Wastewater Samples
 - C. Statistics
 - D. Calculation and Expression of Results
- 7030 COUNTING INSTRUMENTS
 - A. Introduction
 - B. Description and Operation of Instruments
- 7040 FACILITIES
 - A. Counting Room
 - B. Radiochemistry Laboratory
 - C. Laboratory Safety
 - D. Pollution Prevention
 - E. Waste Management
- 7110 GROSS ALPHA AND GROSS BETA RADIOACTIVITY (TOTAL, SUSPENDED, AND DISSOLVED)
 - A. Introduction

- B. Evaporation Method for Gross Alpha-Beta
- C. Coprecipitation Method for Gross Alpha Radioactivity in Drinking Water
- 7120 GAMMA-EMITTING RADIONUCLIDES
 - A. Introduction
 - B. Gamma Spectroscopic Method
- 7500-Cs RADIOACTIVE CESIUM
 - A. Introduction
 - B. Precipitation Method
- 7500-I RADIOACTIVE IODINE
 - A. Introduction
 - B. Precipitation Method
 - C. Ion-Exchange Method
 - D. Distillation Method
- 7500-Ra RADIUM
 - A. Introduction
 - B. Precipitation Method
 - C. Emanation Method
 - D. Sequential Precipitation Method
 - E. Gamma Spectroscopy Method
- 7500-Rn RADON
 - A. Introduction
 - B. Liquid Scintillation Method
- 7500-Sr TOTAL RADIOACTIVE STRONTIUM AND STRONTIUM-90
 - A. Introduction
 - B. Precipitation Method
- 7500-³H TRITIUM
 - A. Introduction
 - B. Liquid Scintillation Spectrometric Method
- 7500-U URANIUM
 - A. Introduction
 - B. Radiochemical Method
 - C. Isotopic Method

Part 8000 TOXICITY

- 8010 INTRODUCTION
 - A. General Discussion
 - B. Terminology
 - C. Basic Requirements
 - D. Conducting Toxicity Tests
 - E. Preparing Organisms for Toxicity Tests
 - F. Toxicity Test Systems, Materials, and Procedures
 - G. Calculating, Analyzing, and Reporting Results of Toxicity Tests
 - H. Interpreting and Applying Results of Toxicity Tests
 - I. Selected Toxicological Literature
- 8020 QUALITY ASSURANCE AND QUALITY CONTROL IN LABORATORY TOXICITY TESTS
 - A. General Discussion

- B. Elements of QA/QC
- 8030 MUTAGENESIS
 - A. Introduction
 - B. *Salmonella* Microsomal Mutagenicity Test
- 8050 BACTERIAL BIOLUMINESCENCE
 - A. Introduction
 - B. Bacterial Bioluminescence Test
- 8070 P450 REPORTER GENE RESPONSE TO DIOXIN-LIKE ORGANIC COMPOUNDS
 - A. Introduction
 - B. The P450 RGS Test
- 8071 COMET/SINGLE-CELL GEL ELECTROPHORESIS ASSAY FOR DETECTION OF DNA DAMAGE
 - A. Introduction
 - B. Comet/Single-Cell Gel Electrophoresis Assay
- 8080 SEDIMENT POREWATER TESTING
 - A. Introduction
 - B. Sediment Collection and Storage
 - C. Extraction of Sediment Pore Water
 - D. Toxicity Testing Procedures
- 8110 ALGAE
- 8111 BIOSTIMULATION (ALGAL PRODUCTIVITY)
 - A. General Principles
 - B. Planning and Evaluating Algal Assays
 - C. Apparatus
 - D. Sample Handling
 - E. Synthetic Algal Culture Medium
 - F. Inoculum
 - G. Test Conditions and Procedures
 - H. Effect of Additions
 - I. Data Analysis and Interpretation
- 8112 PHYTOPLANKTON
 - A. Introduction
 - B. Inoculum
 - C. Test Conditions and Procedures
- 8113 MARINE MACROALGAE
 - A. Introduction
 - B. Selecting and Preparing *Macrocystis pyrifera* Sporophylls
 - C. Toxicity Test Procedures
 - D. Data Evaluation
- 8200 AQUATIC FLOWERING PLANTS
- 8211 DUCKWEED
 - A. Introduction
 - B. Selecting and Preparing Test Organisms
 - C. Toxicity Test Procedure
- 8220 AQUATIC EMERGENT PLANTS
 - A. Introduction
 - B. Selecting and Preparing Test Organisms

- C. Toxicity Test Procedure
- 8310 CILIATED PROTOZOA
 - A. Introduction
 - B. Growth Inhibition Test with Freshwater Ciliate *Colpidium campylum*
 - C. Chemotactic Test with Freshwater Ciliate *Tetrahymena Thermophila*
 - D. Growth Inhibition Test with the Soil Ciliate *Colpoda inflata*
- 8420 ROTIFERS
 - A. Introduction
 - B. Selecting and Preparing Testing Organisms
 - C. Aquatic Toxicity Test Procedures
- 8510 ANNELIDS
 - A. Introduction
 - B. Selecting and Preparing Test Organisms
 - C. Toxicity Test Procedures
 - D. Sediment Test Procedures Using the Marine Polychaete *Neanthes arenaceodentata*
 - E. Sediment Test Procedures using the Marine Polychaete *Polydora cornuta*
 - F. Sediment Test Procedures Using the Freshwater and Marine Oligochaetes *Pristina leidy*, *Tubifex tubifex*, and *Lumbriculus variegatus*
 - G. Data Evaluation
- 8610 MOLLUSKS
 - A. Introduction
 - B. Selecting and Preparing Test Organisms
 - C. Short-Term Test Procedures Using Marine Mollusk Larvae
 - D. Sediment Test Procedures Using Marine Bivalves
 - E. Field Test Procedures Using Freshwater and Marine Bivalves
- 8710 ARTHROPODS
- 8711 *DAPHNIA*
 - A. Introduction
 - B. Selecting and Preparing Test Organisms
 - C. Toxicity Test Procedures
- 8712 *CERIODAPHNIA*
 - A. Introduction
 - B. Selecting and Preparing Test Organisms
 - C. Toxicity Test Procedures
- 8714 MYSIDS
 - A. Introduction
 - B. Selecting and Preparing Test Organisms
 - C. Toxicity Test Procedures
- 8740 DECAPODS
 - A. Introduction
 - B. Selecting and Preparing Test Organisms
 - C. Toxicity Test Procedures
 - D. Data Evaluation
- 8750 AQUATIC INSECTS
 - A. Introduction
 - B. Selecting and Preparing Test Organisms

- C. Toxicity Test Procedures
- D. Data Evaluation
- 8810 ECHINODERM FERTILIZATION AND DEVELOPMENT
 - A. Introduction
 - B. Selecting and Preparing Test Organisms
 - C. Echinoderm Fertilization Test
 - D. Echinoderm Embryo Development Test
- 8910 FISH
 - A. Introduction
 - B. Fish Selection and Culture Procedures
 - C. Toxicity Test Procedures
- 8921 FATHEAD MINNOW
 - A. Introduction
 - B. Culture and Maintenance of Test Organisms
 - C. Toxicity Test Procedures
- 8930 AMPHIBIANS (PROPOSED)
 - A. Introduction
 - B. Culture and Maintenance of Test Organisms
 - C. Toxicity Test Procedures

Part 9000 MICROBIOLOGICAL EXAMINATION

- 9010 INTRODUCTION
- 9020 QUALITY ASSURANCE/QUALITY CONTROL
 - A. Introduction
 - B. Intralaboratory Quality Control Guidelines
 - C. Interlaboratory Quality Control
- 9030 LABORATORY APPARATUS
 - A. Introduction
 - B. Equipment Specifications
- 9040 WASHING AND STERILIZATION
- 9050 PREPARATION OF CULTURE MEDIA
 - A. General Procedures
 - B. Water
 - C. Media Specifications
- 9060 SAMPLES
 - A. Collection
 - B. Preservation and Storage
- 9211 RAPID DETECTION METHODS
 - A. Introduction
 - B. Seven-Hour Fecal Coliform Test
 - C. Special Techniques
- 9212 STRESSED MICROORGANISMS
 - A. Introduction
 - B. Recovery Enhancement
- 9213 RECREATIONAL WATERS
 - A. Introduction

- B. Swimming Pools
 - C. Whirlpools
 - D. Natural Bathing Beaches
 - E. Membrane Filter Technique for *Pseudomonas aeruginosa*
 - F. Multiple-Tube Technique for *Pseudomonas aeruginosa*
- 9215 HETEROTROPHIC PLATE COUNT
- A. Introduction
 - B. Pour Plate Method
 - C. Spread Plate Method
 - D. Membrane Filter Method
 - E. Enzyme Substrate Method
- 9216 DIRECT TOTAL MICROBIAL COUNT
- A. Introduction
 - B. Epifluorescence Microscopic Method Using Acridine Orange
- 9217 ASSIMILABLE ORGANIC CARBON
- A. Introduction
 - B. *Pseudomonas fluorescens* Strain P-17, *Spirillum* Strain NOX Method
- 9218 AEROBIC ENDOSPORES
- A. Introduction
 - B. Membrane Filter Method
- 9221 MULTIPLE-TUBE FERMENTATION TECHNIQUE FOR MEMBERS OF THE COLIFORM GROUP
- A. Introduction
 - B. Standard Total Coliform Fermentation Technique
 - C. Estimation of Bacterial Density
 - D. Presence–Absence (P–A) Coliform Test
 - E. Fecal Coliform Procedure
 - F. *Escherichia coli* Procedure Using Fluorogenic Substrate
 - G. Other *Escherichia coli* Procedures (PROPOSED)
- 9222 MEMBRANE FILTER TECHNIQUE FOR MEMBERS OF THE COLIFORM GROUP
- A. Introduction
 - B. Standard Total Coliform Membrane Filter Procedure
 - C. Delayed-Incubation Total Coliform Procedure
 - D. Thermotolerant (Fecal) Coliform Membrane Filter Procedure
 - E. Delayed-Incubation Thermotolerant (Fecal) Coliform Procedure
 - F. *Klebsiella* Membrane Filter Procedure
 - G. MF Partition Procedures
 - H. Simultaneous Detection of Total Coliform and *E. coli* by Dual-Chromogen Membrane Filter Procedure (PROPOSED)
 - I. Simultaneous Detection of Total Coliform and *E. coli* by Fluorogen/Chromogen Membrane Filter Procedure (PROPOSED)
- 9223 ENZYME SUBSTRATE COLIFORM TEST
- A. Introduction
 - B. Enzyme Substrate Test
- 9224 DETECTION OF COLIPHAGES
- A. Introduction
 - B. Somatic Coliphage Assay

- C. Male-Specific Coliphage Assay Using *Escherichia coli* Famp
 - D. Male-Specific Coliphage Assay Using *Salmonella typhimurium* WG49
 - E. Single-Agar-Layer Method
 - F. Membrane Filter Method
- 9225 DIFFERENTIATION OF COLIFORM BACTERIA
- A. Introduction
 - B. Culture Purification
 - C. Identification
 - D. Media, Reagents, and Procedures
 - E. Reporting Results
- 9230 FECAL STREPTOCOCCUS/STREPTOCOCCUS GROUPS
- A. Introduction
 - B. Multiple-Tube Technique
 - C. Membrane Filter Techniques
 - D. Fluorogenic Substrate Enterococcus Test
- 9240 IRON AND SULFUR BACTERIA
- A. Introduction
 - B. Iron Bacteria
 - C. Sulfur Bacteria
 - D. Enumerating, Enriching, and Isolating Iron and Sulfur Bacteria
 - E. Bacteria Living in Acidic Environments
- 9245 NITRIFYING BACTERIA
- A. Introduction
 - B. Multiple-Tube Method
- 9250 DETECTION OF ACTINOMYCETES
- A. Introduction
 - B. Actinomycete Plate Count
- 9260 DETECTION OF PATHOGENIC BACTERIA
- A. Introduction
 - B. *Salmonella*
 - C. (Reserved)
 - D. (Reserved)
 - E. *Shigella*
 - F. Diarrheagenic *Escherichia coli*
 - G. *Campylobacter*
 - H. *Vibrio*
 - I. *Leptospira*
 - J. *Legionella*
 - K. *Yersinia enterocolitica*
 - L. *Aeromonas*
 - M. *Mycobacterium*
- 9510 DETECTION OF ENTERIC VIRUSES
- A. Introduction
 - B. Virus Concentration from Small Sample Volumes by Adsorption to and Elution from Microporous Filters
 - C. Virus Concentration from Large Sample Volumes by Adsorption to and Elution from

Microporous Filters

- D. Virus Concentration by Aluminum Hydroxide Adsorption-Precipitation
- E. Hydroextraction-Dialysis with Polyethylene Glycol
- F. Recovery of Viruses from Suspended Solids in Water and Wastewater
- G. Assay and Identification of Viruses in Sample Concentrates

9610 DETECTION OF FUNGI

- A. Introduction
- B. Pour Plate Technique
- C. Spread Plate Technique
- D. Membrane Filter Technique
- E. Technique for Yeasts
- F. Zoosporic Fungi
- G. Aquatic Hyphomycetes
- H. Fungi Pathogenic to Humans

9711 PATHOGENIC PROTOZOA

- A. Introduction
- B. Detection of *Giardia* and *Cryptosporidium* in Water
- C. Detection of *Giardia* and *Cryptosporidium* in Wastewater
- D. Infectivity of *Cryptosporidium* in Cell Culture

Part 10000 BIOLOGICAL EXAMINATION

10010 INTRODUCTION

10200 PLANKTON

- A. Introduction
- B. Sample Collection
- C. Concentration Techniques
- D. Preparing Slide Mounts
- E. Microscopes and Calibrations
- F. Phytoplankton Counting Techniques
- G. Zooplankton Counting Techniques
- H. Chlorophyll
- I. Determination of Biomass (Standing Crop)
- J. Metabolic Rate Measurements

10300 PERIPHYTON

- A. Introduction
- B. Sample Collection
- C. Sample Analysis
- D. Primary Productivity
- E. Interpreting and Reporting Results

10400 MACROPHYTES

- A. Introduction
- B. Preliminary Survey
- C. Vegetation Mapping Methods
- D. Population Estimates
- E. Productivity

10500 BENTHIC MACROINVERTEBRATES

- A. Introduction
- B. Sample Collection
- C. Sample Processing and Analysis
- D. Data Evaluation, Presentation, and Conclusions

10600 FISHES

- A. Introduction
- B. Data Acquisition
- C. Sample Preservation
- D. Analysis of Collections
- E. Investigation of Fish Kills

10700 BENTHIC MEIOFAUNA

10750 NEMATOLOGICAL EXAMINATION

- A. Introduction
- B. Collection and Processing Techniques for Nematodes
- C. Illustrated Key to Freshwater Nematodes

10900 IDENTIFICATION OF AQUATIC ORGANISMS

- A. Procedure in Identification
- B. Key to Major Groups of Aquatic Organisms (Plates 1–35)
Acknowledgements
- C. Key for Identification of Freshwater Algae Common in Water Supplies and Polluted Waters
(Plates 1A, 1B, 28–35)
- D. Index to Illustrations
- E. Selected Taxonomic References